

The **GM** Implant



THE GRAND MORSE

Helix GM



NEODENT® GRAND MORSE IMPLANT SYSTEM

GREATNESS IS AN ACHIEVEMENT.

The Neodent® Grand Morse Implant System is the achievement of more than 25 years of experience in implant dentistry, and shared experiences with many clinicians worldwide. Continuing with a unique purpose to always deliver high quality treatment options that changes patients' lives, the Grand Morse Implant System is the Neodent® evolution. Anchor within our philosophy of respecting mechanical and biological principles, this makes it THE implant of choice in dental implant therapy.

 +250
studies

 +25
years
history


 2nd
largest dental implant
company worldwide

The Grand Morse implant system was developed based on the inside out concept, starting from the core of the implant: the prosthetic interface. The result is a solution that combines mechanical strength and versatile prosthetic solutions - from unitary to multiple and from conventional to digital. A complete system that offers several benefits designed to make your work even more efficient.



GRAND RELIABILITY
Stable and strong foundation designed for long term success

GRAND STABILITY
Designed for predictable immediate treatments in all bone types.



GRAND ESTHETICS
Delivers immediate natural esthetics.

GRAND SIMPLICITY
Ease of use at its best.

Grand Morse®

GREATNESS IS AN ACHIEVEMENT



GRAND RELIABILITY

STABLE AND STRONG FOUNDATION DESIGNED FOR LONG TERM SUCCESS

The implant-abutment interface is crucial for a successful long term functional and esthetic result. The Neodent® Grand Morse® connection offers a unique combination based on proven concepts: a platform switching associated with a deep 16° Morse Taper including an internal indexation for a strong and stable connection designed to achieve long-lasting results.



1 Platform Switching

Abutment design with a narrower diameter than the implant coronal area, enabling the platform switching concept⁽⁵⁻⁹⁾.

2 Internal Indexation

Precise abutment positioning, protection against rotation and easy handling.

3 Deep Connection

Allowing a large contact area between the abutment and the implant for an optimal load distribution.

4 16° Morse Taper Connection

Designed to ensure tight fit for an optimal connection sealing.



GRAND SIMPLICITY

EASE OF USE AT ITS BEST

Implant therapy has become an integral part of clinical dentistry, with ever increasing numbers of patients seeking such treatment. The Neodent® Grand Morse® Implant System is smartly engineered providing efficiency and simplicity within the dental treatment network for both surgical to restoratives steps.

ONE PROSTHETIC PLATFORM

All Neodent® Grand Morse® implants feature the unique Grand Morse® connection regardless of the implant diameter.



ONE SCREWDRIVER

The Neo Screwdriver has a star attachment offering reliability and durability compatible with all Neodent® Grand Morse® healing abutments and cover screws and most of the restorative screws.



ONE IMPLANT DRIVER

The Neodent® implant driver allows an easy and reliable implant pick up and placement.



ONE SURGICAL KIT

Intuitive and functional compact surgical kit, that allows the place of Helix GM® implants in all bone types.



GRAND STABILITY

STABLE AND STRONG FOUNDATION DESIGNED FOR LONG TERM SUCCESS

The increasing expectations for shortened treatment duration represent a significant challenge for dental professionals. The Neodent® Grand Morse® system offers a unique implant design featuring the innovative Acqua hydrophilic surface designed to maximize primary stability and predictability in immediate protocols.

HELIX® - OPTIMAL IMPLANT DESIGNED TO ACHIEVE HIGH PRIMARY STABILITY

Helix® Grand Morse® is an innovative hybrid implant design maximizing treatment options and efficiency in all bone types.

Fully tapered body design

- Coronal: 2° - 12°
- Apex: 16°
- » Allowing under-osteotomy




Hybrid contour

- Coronal: Cylindrical
- Apex: Conical
- » For stability with vertical placement flexibility



Active apex

- Soft rounded small tip
- Helical flutes
- » Enabling immediate loading



Dynamic progressive thread design

- Coronal: Trapezoidal > compressing
- Apex: V-Shape > Self-tapping
- » Achieving high primary stability in all bone types



Acqua hydrophilic surface

Designed for high treatment predictability



acqua

Drive®

High primary stability in
challenging bone types.
Bone types III & IV.





GRAND ESTHETICS


DELIVER IMMEDIATE NATURAL ESTHETICS

Nowadays, patients expect both short treatment times and esthetic results. The Neodent®

Grand Morse® restorative portfolio offers flexibility to simplify soft tissue management respecting the biological distances for achieving immediate function and esthetics.



Single-unit screw-retained
prosthesis



Single-unit cement-retained
prosthesis



Overdenture



Multiple-unit screw-retained
prosthesis



Multiple-unit cement-retained
prosthesis



Temporary


NeoPoros

Constant evolution and safety guarantee.

Based on the abrasive sandblasting concept followed by acid etching, the **NeoPoros** surface promotes, by using controlled grain oxides, cavities on the implant surface that then are uniformed with the acid etching technique.

The whole process of obtaining this surface is guaranteed due to automated time, speed, pressure and particle size control.

Several scientific studies continue to be performed so that the **NeoPoros** surface may be always evolving and promoting much more reliability for you.



Controlled roughness on all implant surface. Scanning electron microscopy (A) shows macro (15-30µm) and (B) microtopography (0,3 - 1,3µm).




Image taken by confocal microscopy.
Roughness and Microtopography.
(Sa= 0,3 – 1,3 µm; Sz= 6,0 - 15,5 µm).



Acqua Hydrophilic Surface designed for high treatment predictability.

The Neodent® Acqua hydrophilic surface is the next level of the highly successful S.L.A. type of surface developed to achieve successful outcomes even in challenging situations, such as soft bone or immediate protocols.⁽¹⁻⁴⁾

Hydrophilicity

The hydrophilic surface presents a smaller contact angle when in contact with hydrophilic liquids. This provides greater accessibility of organic fluids to Acqua implant surface.⁽²⁾

Surface comparison

Lab generated images.



NeoPoros surface.



Acqua Hydrophilic
Surface.

Innovative and ease to use

Neodent® Packaging

Neodent® implant packaging has been updated to a concept that provides convenience and safety through all steps of the procedure, from storage to the placement of the implant.

The new packaging aids in identification of both the implant model as well as its diameter and length, regardless of its storage position.



Package instruction of use



After breaking the sterility seal on the blister, hold the primary package (vial) and twist the lid to open it.



To remove the implant from the vial lift the cap up, which has the stand and implant attached to it.



To secure the implant, grip both sides of the implant carrier.



While gripping the implant carrier, remove the lid.



To capture the implant with the contra-angle handpiece attachment, grip the implant carrier while placing the attachment into the implant chamber.





The implant can now be transported to the surgical site.

e-IFU – Electronic Instructions For Use

Neodent® innovates once more, providing an on-line platform designed to provide quick and practical use of its own products instructions: the e-IFU (Instructions For Use) website.

To facilitate access, have the article number, which can be found on the external packaging of the product, in this catalogue or with your local distributor. Once the article number is entered in the website, the professional will have access to relevant information to this product, such as description, indication for use, contraindications, handling, traceability and other features.

Access: ifu.neodent.com.br/en



GROW WITH PEACE OF MIND

Neodent® has developed EasyPack to simplify your daily practice. An all-in-one set that offers everything you need to grow while performing dental implant therapy with confidence, convenience and guidance.

GROW WITH CONFIDENCE
Choose a brand and products you can rely on

GROW WITH CONVENIENCE
The certainty of having everything in one package

GROW WITH GUIDANCE
All workflows in simple steps

THE NEODENT® EASYPACK INCLUDES

- 1 Grand Morse® Helix Implant
- 2 Grand Morse® Cover Screw
- 3 Grand Morse® Healing Abutment
- 4 Grand Morse® Hybrid Implant Analog
- 5 Grand Morse® 3-in-1 Neodent Smart Abutment™ **NEW**





CONVENTIONAL
WORKFLOW



DIGITAL
WORKFLOW

Reliable guided workflow with the 3-in-1 GM Smart Abutment

The combination of the GM Smart Abutment, a unique patented solution combining a closed tray impression coping, a digital scanbody and a temporary abutment in a single piece, with healing components and the analog allows you to choose a restorative path guided for achieving predictable results.



NEODENT® EASYPACK PRODUCT OPTIONS

	Ø 3.5	Ø 3.75	Ø 4.0	Ø 4.3
	Acqua	Acqua	Acqua	Acqua
8.0	138.089	138.113	138.137	138.158
10.0	138.095	138.119	138.143	138.161
11.5	138.101	138.125	138.149	138.164
13.0	138.107	138.131	138.155	138.167
	GM Cover Screw 0 mm	GM Healing Abutment Ø 4.5 X 2.5 mm	GM Hybrid Repositionable Analog* Ø 3.5/3.75 Ø 4.0/4.3 *according to implant diameter	GM Smart Abutment Ø 4.5 X 2.5 mm

Helix GM®

PRODUCT FEATURES:

Implants Description:

- Full dual tapered implant;
- Hybrid contour with a cylindrical coronal part and conical on the apical area;
- Active apex including a soft rounded small tip and helicoidal flutes;
- Dynamic progressive thread design: from compressing trapezoidal threads on the coronal area to self-tapping V-shape threads on the apical part;
- Double threaded implant;
- Grand Morse® connection.

Indications:

- Indicated for all types of bone density and implant immediate placement post extraction.

Drilling features:
























































- Contour drill is required in bone types I and II;
- Final pilot drills are highly recommended in bone types I and II;
- Implant should be positioned 1 or 2 mm below bone level;
- Drilling speed: 800-1200 rpm for bone type I and II;
- Drilling speed: 500-800 rpm for bone type III and IV;
- Implant insertion speed: 30 rpm;
- Maximum torque for implant placement: 60 N.cm.


Available with:




































NeoPoros or 




Drill Sequence
























































																			
	103.170	103.425	103.561	103.578	103.513	103.564	103.579	103.514	103.567	103.580	103.515	103.570	103.581	103.516	103.573	103.582	103.517	103.576	103.577
Ø 3.5																			
Ø 3.75																			
Ø 4.0																			
Ø 4.3																			
Ø 5.0																			


*Optional / Bone types I and II 




































Ø 3.5																			
Ø 3.75							*												
Ø 4.0									*										
Ø 4.3													*						
Ø 5.0																	*		
Ø 6.0																			
Ø 7.0																			


*Optional / Bone types III and IV 

Drill Sequence with Neodent® Control System









																			
	103.170	103.492	103.493	103.500	103.513	103.494	103.501	103.514	103.495	103.502	103.515	103.496	103.503	103.516	103.497	103.504	103.517	103.498	103.499
Ø 3.5																			
Ø 3.75																			
Ø 4.0																			
Ø 4.3																			
Ø 5.0																			

*Optional / Bone types I and II 


Ø 3.5																			
Ø 3.75							*												
Ø 4.0									*										
Ø 4.3													*						
Ø 5.0																	*		
Ø 6.0																			
Ø 7.0																			

*Optional / Bone types III and IV 


Helix GM® Implants

Ø 3.5	Acqua	NeoPoros	Ø 3.75	Acqua	NeoPoros	Ø 4.0	Acqua	NeoPoros	Ø 4.3	Acqua	NeoPoros
	8.0 140.943	109.943		8.0 140.976	109.976		8.0 140.982	109.982		8.0 140.948	109.948
	10.0 140.944	109.944		10.0 140.977	109.977		10.0 140.983	109.983		10.0 140.949	109.949
	11.5 140.945	109.945		11.5 140.978	109.978		11.5 140.984	109.984		11.5 140.950	109.950
	13.0 140.946	109.946		13.0 140.979	109.979		13.0 140.985	109.985		13.0 140.951	109.951
	16.0 140.947	109.947		16.0 140.980	109.980		16.0 140.986	109.986		16.0 140.952	109.952
	18.0 140.988	109.988		18.0 140.981	109.981		18.0 140.987	109.987		18.0 140.989	109.989
Ø 5.0	Acqua	NeoPoros	Ø 6.0	Acqua	NeoPoros	Ø 7.0	Acqua	NeoPoros	GM Cover Screw		
	8.0 140.953	109.953		8.0 140.1009	109.1009		8.0 140.1059	109.1059			
	10.0 140.954	109.954		10.0 140.1010	109.1010		10.0 140.1060	109.1060	0 mm 117.021		
	11.5 140.955	109.955		11.5 140.1011	109.1011		11.5 140.1061	109.1061	2 mm 117.022		
	13.0 140.956	109.956		13.0 140.1012	109.1012		13.0 140.1062	109.1062	:: Use the manual Neo Screwdriver (104.060);		
	16.0 140.957	109.957							:: Do not exceed the insertion torque of 10 N.cm.		
	18.0 140.990	109.990									

GM Healing Abutment

	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm
	Ø 3.3 106.207	106.208	106.209	106.210	106.211	106.212
	Ø 4.5 106.213	106.214	106.215	106.216	106.217	106.218
:: Use the manual Neo Screwdriver (104.060);						
:: Do not exceed the insertion torque of 10 N.cm.						

GM Customizable Healing Abutment

	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm	6.5 mm
	Ø 5.5 106.223	106.224	106.225	106.226	106.227	
	Ø 7.0 106.228	106.229	106.230	106.231	106.232	
:: Use the manual Neo Screwdriver (104.060);						
:: Do not exceed the insertion torque of 10 N.cm.						

Drive GM®

PRODUCT FEATURES:

Implants Description:

- Tapered implant;
- Square shape threads;
- Double threaded implant;
- Reverse cutting chambers distributed across the implant body;
- Rounded apex with a sharp edge;
- Grand Morse® connection.

Indications:

- Indicated for bone types III and IV and implant immediate placement post-extraction;

Drilling features:

- Final pilot drill is optional in bone types III and IV;
- Implant should be positioned 1 or 2 mm below bone level;
- Drilling speed: 500-800 rpm;
- Implant insertion speed: 30 rpm;
- Maximum torque for implant placement: 60 N.cm.



Available with:



Drill Sequence

	Initial	Ø 2.0	Ø 3.5	Ø 3.5	Ø 4.3	Ø 4.3	Ø 5.0	Ø 5.0
	103.170	103.425	103.561	103.513	103.570	103.516	103.573	103.517
Ø 3.5 mm								
Ø 4.3 mm								
Ø 5.0 mm								

*Optional / Bone types III and IV

Drive GM® Implants

		8.0 mm	10.0 mm	11.5 mm	13.0 mm	16.0 mm	18.0 mm
Ø 3.5							
Acqua	140.958	140.959	140.960	140.961	140.962	140.963	
Ø 4.3							
Acqua	140.964	140.965	140.966	140.967	140.968	140.969	
Ø 5.0							
Acqua	140.970	140.971	140.972	140.973	140.974	140.975	

GM Healing Abutment

	Profile	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm
Ø 3.3	106.207	106.208	106.209	106.210	106.211	106.212	
Ø 4.5	106.213	106.214	106.215	106.216	106.217	106.218	

:: Use the manual Neo Screwdriver (104.060);
:: Do not exceed the insertion torque of 10 N.cm.

GM Cover Screw


		0 mm	2 mm
	117.021	117.022	

:: Use the manual Neo Screwdriver (104.060);
:: Do not exceed the insertion torque of 10 N.cm.

GM Customizable Healing Abutments

	Profile	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm	6.5 mm
Ø 5.5	106.223	106.224	106.225	106.226	106.227		
Ø 7.0		106.228	106.229	106.230	106.231	106.232	

GM Abutment



Recommended for posterior region.

Consider in addition 1.5 - 2.0 mm
for the restorative material;

Minimum interocclusal space of 4.9
mm from the mucosa level;

Exact;

Unlocking feature.



GM Anatomic Abutment



Recommended for anterior region.

Gingiva color for esthetic outcomes;

Click retention for provisional copings;

Exact;

Unlocking feature.



Installation Sequence

0.8 mm	1.5 mm	2.5 mm	GM Exact Abutment
115.237	115.238	115.239	
3.5 mm	4.5 mm	5.5 mm	
115.240	115.241	115.242	



Drivers




Accessories



*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

Installation Sequence


In Mouth




Drivers




In Lab



GM Mini Conical Abutment



Multiple-unit screw-retained prosthesis




Ø 4.8 mm


Consider in addition 1.5 - 2.0 mm for the restorative material;
Minimum interocclusal space of 4.5 mm from the mucosa level for straight abutments.




GM Micro Abutment



Single-unit screw-retained prosthesis



Multiple-unit screw-retained prosthesis




Ø 3.5 mm

Consider in addition 1.5 - 2.0 mm for the restorative material;
Minimum interocclusal space of 3.5 mm from the mucosa level.



Recommended for limited spaces and narrow inter-dental spaces.

Installation Sequence



Drivers

1 Hexagonal Prosthetic Driver + Torque Wrench

2 Neo Screwdriver Torque Connection + Torque Wrench

3 Neo Screwdriver Torque Connection + Manual Screwdriver Torque


Accessories

Mini Conical Abutment Polishing Protector 123.008

Replacement Coping Screw 116.269 Titanium 116.270 Neotorque*

*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

Installation Sequence



Drivers

1 Hexagonal Prosthetic Driver + Torque Wrench

2 Neo Screwdriver Torque Connection + Torque Wrench

3 Neo Screwdriver Torque Connection + Manual Screwdriver Torque

Accessories

Micro Abutment Polishing Protector 123.015 Bridge


Replacement Coping Screw 116.269 Titanium 116.270 Neotorque*

*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

GM Universal Abutment

Single-unit cement-retained prosthesis


Ø 3.3/4.5 mm




GM Titanium Base

Single-unit screw-retained prosthesis


Single-unit cement-retained prosthesis

Ø 3.5/4.5/5.5/6.5 mm


With removable screw.



Installation Sequence



Installation Sequence




Drivers



Drivers




Accessories



*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

GM Pro Peek Abutment



Single-unit cement-retained temporary prosthesis



Ø 4.5/ 6.0 mm

Biocompatible Peek of easy customization.

Consider in addition 1.5 - 2.0 mm for the restorative material;


Interocclusal height of 9.2 mm (can be customized up to 5.0 mm);

Exact;

Unlocking feature.




Installation Sequence



GM Titanium Base for Bridge



Multiple-unit screw-retained prosthesis



Multiple-unit cement-retained prosthesis




Ø 3.5/4.5/ 5.5 mm

With removable screw.


Cementable area:
4.0 mm for Ø 3.5
4.5 mm for Ø 4.5
and Ø 5.5.



Installation Sequence




Drivers



Drivers



Accessories



*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

Titanium Base C for GM


Single-unit screw-retained prosthesis

Single-unit cement-retained prosthesis

Ø 4.65 mm

With removable screw.

Installation Sequence



Workflow

Step 1
Gingiva height selection and ordering.

Select the Titanium Base C for GM Exact gingival height.

Step 2
Intra-oral scanning.

Insert the Titanium Base C for GM Exact in the Neodent® implant.

Step 3
Design and milling.

Select in the CAD software the comparable third-party Ti-base and perform the digital design.

Step 4
Finalization and fixation.

Check the fit of milled restoration in the patient's mouth and adapt it, if needed.
Cement the restoration on the Titanium Base C for GM Exact and insert it into the patient's mouth.

Order the Titanium Base C for GM Exact.
Please note that the scanbody has to be purchased directly from equipment manufacturer.

Insert scanbody on the Titanium Base C for GM Exact.

Mill the digital design.

CEREC digital library compatibility

Library	Sirona's Products				Compatible with implant System	
Ti-base	Scanbody	REF Scanbody Omnicam	REF Scanbody Bluecam / Ineos	Grinding block	Implant manufacturer	Implant system
NBB 3.4 L						
NB A 4.5 L						
SSO 3.5 L						
S BL 3.3 L						
S BL 4.1 L						
BO 3.4 L						

Drivers

1 Neo Screwdriver Torque Connection

+

Torque Wrench

Accessories

Replacement Sterile Screws

116.286 Titanium

116.285 Neotorque*

*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

GM Titanium Base Angled Solution (AS)


Single-unit screw-retained prosthesis

Single-unit cement-retained prosthesis

Ø 4.0/4.5/5.5 mm

With removable screw.

Installation Sequence



0.8 mm 1.5 mm 2.5 mm

Ø 4.0 135.327 135.328 135.329

Ø 4.5 135.333 135.334 135.335

Ø 5.5 135.339 135.340 135.341

GM Titanium Base Angled Solution (AS) 4mm

or

GM Titanium Base Angled Solution (AS) 6mm

0.8 mm 1.5 mm 2.5 mm

Ø 4.0 135.330 135.331 135.332

Ø 4.5 135.336 135.337 135.338

Ø 5.5 135.342 135.343 135.344

Drivers

1 Angled Solution Screwdriver for Torque Wrench

105.150 Short 105.151 Regular 105.152 Long

+

Torque Wrench

or

1 Angled Solution Screwdriver for Contra-angle

105.147 Short 105.148 Regular 105.149 Long

+

Contra-angle

Accessories

2 Neo Screwdriver Torque Connection


+


Manual Screwdriver Torque


Replacement Sterile Screw

116.288 Screw for GM Titanium Base AS


GM Titanium Block for MEDENTiKA Holder

Single-unit screw-retained prosthesis

Single-unit cement-retained prosthesis


Multiple-unit cement-retained prosthesis


Ø 11.5/
15.8 mm




Screw sold separately.

GM Titanium Block for AG Holder

Single-unit screw-retained prosthesis

Single-unit cement-retained prosthesis

Multiple-unit cement-retained prosthesis

Ø 12.0 mm




Screw sold separately.

Installation Sequence

Complete Digital Workflow



Semi Digital Workflow




Installation Sequence

Complete Digital Workflow



Semi Digital Workflow



Drivers

1 Neo Screwdriver Torque Connection



+



2 Neo Screwdriver Torque Connection

Manual Screwdriver Torque

Accessories



Sterile Screws sold separately

116.286 Titanium

116.285 Neotorque*

*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

Drivers

1 Neo Screwdriver Torque Connection




+



2 Neo Screwdriver Torque Connection

Manual Screwdriver Torque

Accessories




Sterile Screws sold separately


116.286 Titanium


116.285 Neotorque*

*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

GM CoCr Abutment

Single-unit screw-retained prosthesis

Single-unit cement-retained prosthesis

Ø 4.1/4.5/5.0 mm

For implants placed at bone level.

Consider in addition 1.5 - 2.0 mm for the restorative material;
Interocclusal height of 12 mm (can be customized up to 5.0 mm);

Exact.



GM Temporary Abutment

Single-unit screw-retained temporary prosthesis

Multiple-unit screw-retained temporary prosthesis


Ø 3.5/4.5 mm

Customizable area made of titanium.
A minimum height of 4 mm of the customizable area must be kept.
With retentive grooves for acrylic material and allows customization.

Consider in addition 1.5 - 2.0 mm for the restorative material;
Channels of customizations;
Interocclusal height of 10 mm (can be customized up to 4.0 mm);
Exact.



Installation Sequence



Drivers

Neo Screwdriver Torque Connection

+

Torque Wrench

Neo Screwdriver Torque Connection

+

Manual Screwdriver Torque


Accessories

Replacement Sterile Screws

116.283 Titanium
116.282 Neotorque*

*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

Installation Sequence



Drivers

Neo Screwdriver Torque Connection

+

Torque Wrench


Accessories

Replacement Sterile Screws

116.286 Titanium
116.285 Neotorque*

*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

GM Novaloc



Angled version with removable screw.




Grand Morse® Surgical Kit

Autoclavable polymer case.
To order the pre-mounted version of the kit, with its complete composition, use code [110.302](#).



Installation Sequence



Articles

- 110.288 GM Surgical Kit Case
- 103.162 Twist Drill 2.0 Plus
- 103.213 Pilot Drill 2.0/3.0 Plus
- 103.164 Twist Drill 3.0 Plus
- 103.166 Twist Drill 3.3 Plus
- 103.167 Twist Drill 3.8 Plus
- 103.168 Twist Drill 4.3 Plus
- 103.163 Twist Drill 2.8 Plus
- 103.170 Initial Drill Plus
- 103.414 Pilot Drill GM 2.8/3.5
- 103.415 Pilot Drill GM 3.0/3.75
- 103.416 Pilot Drill GM 3.3/4.0
- 103.417 Pilot Drill GM 4.3
- 103.418 Pilot Drill GM 4.3/5.0

- 103.419 Tapered Contour Drill 3.5
- 103.420 Tapered Contour Drill 3.75
- 103.421 Tapered Contour Drill 4.0
- 103.422 Tapered Contour Drill 4.3
- 103.423 Tapered Contour Drill 5.0
- 103.425 Tapered Drill 2.0
- 103.399 Tapered Drill 3.5
- 103.402 Tapered Drill 3.75
- 103.405 Tapered Drill 4.0
- 103.408 Tapered Drill 4.3
- 103.411 Tapered Drill 5.0
- 103.427 Tapered Drill 6.0
- 105.131 GM Implant Driver - Contra-Angle
- 104.060 Neo Screwdriver (Medium)

- 105.130 GM Implant Driver - Torque Wrench (Long)
- 104.028 Manual Implant Driver - Contra-Angle
- 105.129 GM Implant Driver - Torque Wrench (Short)
- 128.019 Direction Indicator 2.8/3.5
- 128.020 Direction Indicator 3.0/3.75
- 128.021 Direction Indicator 3.3/4.0
- 128.022 Direction Indicator 3.6/4.3
- 128.023 Direction Indicator 4.3/5.0
- 128.028 Height Mesurer GM
- 129.004 Depth Probe
- 129.001 Titanium Tweezers
- 104.050 Torque Wrench
- 103.426 Drill Extension

Note: Items that compose Neodent® Kits are sold separately.

Helix GM® Compact Surgical Kit

Autoclavable polymer case.
The Kit allows the installation of Helix GM® Implants in all bone types.
To order the pre-mounted version of the kit, with its complete composition, use code [110.303](#).



Articles

- 110.297 Helix GM® Compact Surgical Kit Case
- 103.170 Initial Drill
- 103.425 Tapered Drill 2.0
- 103.399 Tapered Drill 3.5
- 103.402 Tapered Drill 3.75
- 103.405 Tapered Drill 4.0
- 103.408 Tapered Drill 4.3
- 103.411 Tapered Drill 5.0
- 103.427 Tapered Drill 6.0
- 103.487 Tapered Drill 7.0 (Short)*
- 104.060 Neo Manual Screwdriver (Medium)
- 104.028 Manual Implant Driver - Contra-angle


- 103.426 Drill Extension
- 103.419 Tapered Contour Drill 3.5
- 103.420 Tapered Contour Drill 3.75
- 103.421 Tapered Contour Drill 4.0
- 103.422 Tapered Contour Drill 4.3
- 103.423 Tapered Contour Drill 5.0
- 105.131 GM Implant Driver - Contra-angle GM
- 105.130 Implant Driver - Torque Wrench (Long)
- 105.129 GM Implant Driver - Torque Wrench (Short)
- 103.414 GM Pilot Drill 2.8/3.5
- 103.415 GM Pilot Drill 3.0/3.75
- 103.416 GM Pilot Drill 3.3/4.0

- 103.417 GM Pilot Drill 4.3
- 103.418 GM Pilot Drill 4.3/5.0
- 128.028 GM Height Mesurer
- 128.030 Angle Mesurer for Drill 2.0 17°
- 128.031 Angle Mesurer for Drill 2.0 30°
- 128.019 Direction Indicator 2.8/3.5
- 128.020 Direction Indicator 3.0/3.75
- 128.021 Direction Indicator 3.3/4.0
- 128.022 Direction Indicator 3.6/4.3
- 128.023 Direction Indicator 4.3/5.0
- 129.004 Depth Probe
- 104.050 Torque Wrench

Note: Items that compose Neodent® Kits are sold separately.

*Tapered Drill 7.0 is not included in the pre-mounted kit composition (110.303).


Drivers






Accessories




Measurements GM
Mini Conical Abutment



Measurements GM
Anatomic Abutment



Measurements GM
Universal Abutment



Neodent controlsystem

User friendly kit retentive system

The Neodent® Control Drill Stop Kit includes an innovative retentive system.



TRUST YOURSELF

The surgical procedure for implant placement can be perceived as complex, especially when performed in the posterior regions with limited visibility, or in proximity with anatomical structures such as nerve canals. The Neodent® Control System brings confidence and efficiency building trust during the surgical procedure.

Protect anatomical structures

The placement of implants requires accuracy, and the Neodent® Control System has been designed to reduce the risk against overdrilling and protecting anatomical structures such as nerves, the sinus or adjacent roots by securing the final depth.

Master limited visibility

The Neodent® Control System helps to provide confidence during situations with reduced visibility due to adjacent teeth, limited mouth opening, blood, saliva, making it difficult to read the lines on a twisting drill by reaching the planned depth.



Intuitive solution

The Neodent® Control System is a color coded solution facilitating the identification of the drill sequence, the diameter and length of the implant and the combination of drill stop and drill.

Secure drill stop locking system

The Neodent® Control Drill Stop features a modern drill locking system enabling an easy and secure engaging into the drill, offering a peace-of-mind surgical experience.

Multiple use solution

The Neodent® Control Drill Stops are made of titanium for professional cleaning and autoclaving allowing multiple use.



A convenient and time-saving pick and drop mechanism during the surgical procedure.

Neodent® Color Code overview



Color code according to implant length



Compatible portfolio of Helix GM® Implants



Length	Diameter						
	3.5	3.75	4.0	4.3	5.0	6.0	7.0
8	✓	✓	✓	✓	✓	✓	✓
10	✓	✓	✓	✓	✓	✓	✓
11.5	✓	✓	✓	✓	✓	✓	✓
13	✓	✓	✓	✓	✓	✓	✓

Helix GM® Compact Kit Control Stop Drills

Autoclavable polymer case.
The Kit allows the installation of Helix GM® Implants in all bone types, using the Neodent® Control Stop Drills.
To order the pre-mounted version of the kit, with its complete composition, use code [110.308](#).



Articles

- 110.297

Helix GM® Compact Surgical Kit Case
- 103.170

Initial Drill
- 103.492

Tapered Control Stop Drill 2.0
- 103.493

Tapered Control Stop Drill 3.5
- 103.494

Tapered Control Stop Drill 3.75
- 103.495

Tapered Control Stop Drill 4.0
- 103.496

Tapered Control Stop Drill 4.3
- 103.497

Tapered Control Stop Drill 5.0
- 103.498

Tapered Control Stop Drill 6.0 (Short)
- 103.499

Tapered Control Stop Drill 7.0 (Short)*
- 104.060

Neo Manual Screwdriver (Medium)
- 104.028

Manual Implant Driver - Contra-angle

103.426

Drill Extension

103.500

Tapered Control Stop Drill 3.5+

103.501

Tapered Control Stop Drill 3.75+

103.502

Tapered Control Stop Drill 4.0+

103.503

Tapered Control Stop Drill 4.3+

103.504

Tapered Control Stop Drill 5.0+

105.131

GM Implant Driver - Contra-angle GM

105.130

Implant Driver - Torque Wrench (Long)

105.129

GM Implant Driver - Torque Wrench (Short)

103.513

Pilot Drill 3.5

103.514

Pilot Drill 3.75

103.515

Pilot Drill 4.0

- 103.516

Pilot Drill 4.3
- 103.517

Pilot Drill 5.0
- 128.028

GM Height Measurer
- 128.030

Angle Measurer for Drill 2.0 17°
- 128.031

Angle Measurer for Drill 2.0 30°
- 128.019

Direction Indicator 2.8/3.5
- 128.020

Direction Indicator 3.0/3.75
- 128.021

Direction Indicator 3.3/4.0
- 128.022

Direction Indicator 3.6/4.3
- 128.023

Direction Indicator 4.3/5.0
- 129.004


Depth Probe
- 104.050

Torque Wrench

Note: Items that compose Neodent® Kits are sold separately.
*Tapered Control Stop Drill 7.0 is not included in the pre-mounted kit composition (110.308).

Grand Morse® Prosthetic Kit

Autoclavable polymer case.
To order the pre-mounted version of the kit, with its complete composition, use code [110.304](#).



Articles

- 110.294

GM Prosthetic Kit Case
- 105.146

Neo Screwdriver Torque Connection - Contra-angle (Extra-short)
- 105.135

Neo Screwdriver Torque Connection - Contra-angle (Short)
- 105.136

Neo Screwdriver Torque Connection - Contra-angle (Medium)
- 105.138

Hexagonal Prosthetic Driver - Contra-angle
- 105.137

Hexagonal Prosthetic Driver - Torque Wrench
- 105.133

Neo Screwdriver Torque Connection (Short) - Torque Wrench
- 105.132

Neo Screwdriver Torque Connection (Medium) - Torque Wrench
- 105.134

Neo Screwdriver Torque Connection (Long) - Torque Wrench
- 104.005

Manual Screwdriver Torque
- 128.028

GM Height Measurer
- 104.050

Torque Wrench

Note: Items that compose Neodent® Kits are sold separately.

Control Drill Stop Kit

Autoclavable polymer case.
The Kit allows the sterilization and engagement of Neodent® Control Drill Stops on the drills.
To order the pre-mounted version of the kit, with its complete composition, use code [110.306](#).



Articles

- 110.307

Control Drill Stop Kit Case
- 125.144

8.0 Control Drill Stop D2.0
- 125.145

10.0 Control Drill Stop D2.0
- 125.146

11.5 Control Drill Stop D2.0
- 125.147

13.0 Control Drill Stop D2.0
- 125.148

8.0 Control Drill Stop D3.5
- 125.149

10.0 Control Drill Stop D3.5
- 125.150

11.5 Control Drill Stop D3.5
- 125.151

13.0 Control Drill Stop D3.5
- 125.152

8.0 Control Drill Stop D3.75/4.0
- 125.153

10.0 Control Drill Stop D3.75/4.0
- 125.154

11.5 Control Drill Stop D3.75/4.0

125.155

13.0 Control Drill Stop D3.75/4.0

125.156

8.0 Control Drill Stop D4.3/5.0

125.157

10.0 Control Drill Stop D4.3/5.0

125.158

11.5 Control Drill Stop D4.3/5.0

125.159

13.0 Control Drill Stop D4.3/5.0

125.160

8.0 Control Drill Stop D6.0/7.0

125.161

10.0 Control Drill Stop D6.0/7.0

125.162

11.5 Control Drill Stop D6.0/7.0

125.163

13.0 Control Drill Stop D6.0/7.0

Note: Items that compose Neodent® Kits are sold separately.

Grand Morse® Try-In Kit

Autoclavable polymer case.
To order the pre-mounted version of the kit, with its complete composition, use code [110.305](#).



Articles

- 110.295

GM Try-In Kit Case
- 114.772

GM Abutment Try-In 3.3X6X0.8
- 114.773

GM Abutment Try-In 3.3X6X1.5
- 114.774

GM Abutment Try-In 3.3X6X2.5
- 114.775

GM Abutment Try-In 3.3X6X3.5
- 114.776

GM Abutment Try-In 3.3X6X4.5
- 114.777

GM Abutment Try-In 3.3X6X5.5
- 114.778

GM Abutment Try-In 4.5X6X0.8
- 114.779

GM Abutment Try-In 4.5X6X1.5
- 114.780

GM Abutment Try-In 4.5X6X2.5
- 114.781

GM Abutment Try-In 4.5X6X3.5

114.782

GM Abutment Try-In 4.5X6X4.5

114.783

GM Abutment Try-In 4.5X6X5.5

114.784

GM Abutment Try-In 17° 3.3X6X1.5

114.785

GM Abutment Try-In 17° 3.3X6X2.5

114.786

GM Abutment Try-In 17° 3.3X6X3.5

114.787

GM Abutment Try-In 17° 4.5X6X1.5

114.788

GM Abutment Try-In 17° 4.5X6X2.5

114.789

GM Abutment Try-In 17° 4.5X6X3.5

114.790

GM Abutment Try-In 30° 3.3X6X1.5

114.791

GM Abutment Try-In 30° 3.3X6X2.5

114.792

GM Abutment Try-In 30° 3.3X6X3.5

114.793

GM Abutment Try-In 30° 4.5X6X1.5

114.794

GM Abutment Try-In 30° 4.5X6X2.5

114.795

GM Abutment Try-In 30° 4.5X6X3.5

114.796

GM Anatomic Abutment Try-In 1.5

114.797

GM Anatomic Abutment Try-In 2.5

114.798

GM Anatomic Abutment Try-In 3.5

114.799

GM Lateral Anatomic Abutment Try-In 1.5

114.800

GM Lateral Anatomic Abutment Try-In 2.5

114.801

GM Lateral Anatomic Abutment Try-In 3.5

104.058

Neo Manual Screwdriver (Short)

128.028

GM Height Measurer

Note: Items that compose Neodent® Kits are sold separately.

Grand Morse® Instruments



Initial Drill

- :: Available in surgical steel;
- :: 2.0mm diameter.

103.170

Tapered Drills

- :: Available in surgical steel;
- :: Drill sequence for Helix GM® and Drive GM® Implants;
- :: With a color code according to the drill diameter.



	Short 31 mm	Regular 35 mm	Long 43 mm
Ø 2.0	103.559	103.425	103.560
Ø 3.5	103.562	103.561	103.563
Ø 3.75	103.565	103.564	103.566
Ø 4.0	103.568	103.567	103.569
Ø 4.3	103.571	103.570	103.572
Ø 5.0	103.574	103.573	103.575
Ø 6.0	103.576		
Ø 7.0	103.577		

Tapered+ Drills

- :: For preparing the implant bed in bone types I and II for Helix GM® Implants;
- :: With a color code according to the drill diameter and 2 stripes of color for identification.



Ø 3.5+	103.578
Ø 3.75+	103.579
Ø 4.0+	103.580
Ø 4.3+	103.581
Ø 5.0+	103.582

Pilot Drills

- :: Available in surgical steel;
- :: Increasing the surgical alveolus diameter ridge, easing the penetration of the next drill or the implant.



Ø 2/3	103.213		
Ø 3.5	103.513	Ø 5.0	103.517
Ø 3.75	103.514	Ø 3.8/4.3	103.214
Ø 4.0	103.515	Ø 4.3/5.3	103.215
Ø 4.3	103.516	Ø 5.3/6	103.221



Tapered Control Stop Drills

- :: Available in surgical steel;
- :: Drill sequence for Helix GM® Implants;
- :: Attachment to engage drill stops;
- :: With a color code according to the drill diameter.

Ø 2.0	103.492	Ø 4.3	103.496
Ø 3.5	103.493	Ø 5.0	103.497
Ø 3.75	103.494	Ø 6.0	103.498
Ø 4.0	103.495	Ø 7.0	103.499



Tapered+ Control Stop Drills

- :: Available in surgical steel;
- :: For preparing the implant bed in bone types I and II for Helix GM® Implants;
- :: Attachment to engage drill stops;
- :: With a color code according to the drill diameter and 2 stripes of color for identification.

Ø 3.5+	103.500	Ø 4.3+	103.503
Ø 3.75+	103.501	Ø 5.0+	103.504
Ø 4.0+	103.502		

Control Drill Stops

- :: Available in titanium;
- :: To be used in association with the Control Stop Drills;
- :: Physical control for drilling depth.




	8 mm	10 mm	11.5 mm	13 mm
Ø 2.0	125.144	125.145	125.146	125.147
Ø 3.5	125.148	125.149	125.150	125.151
Ø 3.75/4.0	125.152	125.153	125.154	125.155
Ø 4.3/5.0	125.156	125.157	125.158	125.159
Ø 6.0/7.0	125.160	125.161	125.162	125.163



Direction Indicators

- :: Available in titanium;
- :: Instrument to guide the implant position;
- :: Diameter of central band corresponds to GM Implant diameter;
- :: Smaller side to be used after Ø2.0mm drill;
- :: Larger side to be used after the last drill before implant installation.



2.8/3.5	128.019	3.6/4.3	128.022
3.0/3.75	128.020	4.3/5.0	128.023
3.3/4.0	128.021		



Drill Extension

- :: Available in surgical steel;
- :: Fit the drill directly into the Drill Extension.

103.426



GM Height Measurer

- :: Available in titanium;
- :: For selecting GM prosthetic abutments;
- :: Marks corresponding to transmucosa heights.
- :: Can be used as X-Ray Positioner.

128.028

GM Implant Driver - Contra-Angle

- :: To capture the implant directly from the packaging;
- :: To place GM Implants with contra-angle, or attached to a manual driver for contra-angle connections (104.028) for hand placement;
- :: With six dimples to indicate the hex index face position;
- :: The laser marks indicate the depth of implant placement, bone level, 1 and 2mm infra-bone and last marking (3mm) biological space;
- :: Maximum torque 35 N.cm.

105.131

GM Implant Driver - Torque Wrench

- :: To place GM Implants with the Torque Wrench (104.050);
- :: With six marks to indicate the hex index face position;
- :: The laser marks indicate the depth of implant placement, bone level, 1 and 2mm infra-bone and last marking (3mm) biological space;
- :: Maximum torque: 60 N.cm..

Short 22 mm	Long 30 mm
105.129	105.130

Neo Screwdriver Torque Connection - Torque Wrench

- :: Available in surgical steel;
- :: Yellow color for line identification.

Short 16.5 mm	Medium 22 mm	Long 32 mm
105.133	105.132	105.157

Neo Manual Screwdriver

- :: Available in surgical steel;
- :: Yellow color for line identification

Short 21 mm	Medium 25 mm	Long 37 mm
104.058	104.060	104.072

Neo Screwdriver Torque Connection - Contra-angle

- :: Available in surgical steel;
- :: Yellow color for line identification;
- :: Extra Short Neo Screwdriver Torque Connection - Contra-angle (105.146) recommended for Impression Copings, Cover Screws and Healing Abutments.

Extra Short 16.5 mm	Short 24 mm	Long 31 mm
105.146	105.135	105.160

Hexagonal Prosthetic Driver

- :: Available in surgical steel;
- :: To install and apply torque over straight GM Mini Conical Abutments and GM Micro Abutments;

Contra-angle	Torque Wrench
105.138	105.137

Angled Solution Screwdriver for Torque Wrench

- :: To place GM Titanium Bases for Angled Solution with torque wrench;
- :: Maximum torque of 20 N.cm.

Short 16.5 mm	Medium 22.5 mm	Long 28.5 mm
105.150	105.151	105.152

Angled Solution Screwdriver for Contra-angle

- :: To place GM Titanium Bases for Angled Solution with contra-angle;
- :: Maximum torque of 20 N.cm.

Short 20 mm	Medium 26 mm	Long 32 mm
105.147	105.148	105.149

GM Bone Profile Drill with Guide

- :: Available in surgical steel;
- :: Used in the surgical second step;
- :: Conforms the bone around the implant platform, preparing the emergence profile to be suitable to prosthetic components.

103.424

Angle Measurer for Drill 2.0

- :: Available in titanium;
- :: Angles: 17° and 30°;
- :: To select and plan the abutments angulation during surgical procedures;
- :: Suggested use: after Twist Drill 2.0.

17°	30°
128.030	128.031

GM Angle Measurer

- :: Available in titanium;
- :: Angles: 17° and 30°;
- :: To a more accurate selection and planning of the abutments angulation during the prosthetic phase.

17°	30°
128.032	128.033

Control Stop Kit Holder

- :: Available in polymer;
- :: Replacement piecel;
- :: To keep the stops organized and to engage and remove them from the drills.

110.310

Manual Implant Drivers

- :: Available in surgical steel;
- :: For Contra-angle connections: connected to GM Implant Driver, it becomes a manual driver for implant placement.
- :: For Torque Wrench connections: connected to screwdrivers, it provides manual torque.

Contra-angle Connections	Torque Wrench Connections
104.028	104.005

Remover for Abutments with internal threads

- :: Available in surgical steel;
- :: To remove abutments with internal threads from the implants, after removal of the screws;
- :: Compatible with abutments with Neo removable Screws

130.118	Long 130.114
---------	-----------------

Remover for Neo Screws

- :: Available in surgical steel;
- :: Compatible with Neo remvoable screws for abutments

130.119	Long 130.115
---------	-----------------

Torque Wrench

- :: Available in surgical steel;
- :: Fitting for square connections;
- :: Collapsible Wrench that allows for proper assembly cleaning.

104.050

Removal Sets for Abutments with internal threads and Neo Screws

- :: Available in surgical steel;
- :: To remove Neo Removable Screws and abutments with internal threads from the implants, after removal of the screws;
- :: Compatible with abutments with Neo removable Screws

130.117	Long 130.116
---------	-----------------

SIMPLICITY AT ONE HAND

Neodent® is designed to offer straightforward guided surgery techniques enabling predictable surgical results, efficient treatment protocols and patient treatment acceptance.

- ✓ **STRAIGHTFORWARD GUIDED SURGERY TECHNIQUE**
Surgical convenience with one-hand procedures
- ✓ **EFFICIENT TREATMENT PROTOCOLS**
Intuitive and simple technique
- ✓ **PREDICTABLE SURGICAL RESULTS**
Confidence for accurate implant positioning
- ✓ **PATIENT TREATMENT ACCEPTANCE**
Communication building trust and patient engagement




NEODENT® EASYGUIDE ENABLES ONE-HAND PROCEDURES WITH NO DRILL HANDLES

- Simple technique
- Reduced number of instruments
- Surgeries can be performed without assistance

ONE DRILL DESIGN

The unique geometry of the Neodent® **EasyGuide** tapered drills is indicated for all bone types and dismisses the need for additional drill types or taps, simplifying the drilling sequence.



FULLY GUIDED BED PREPARATION

- Intimate contact between drill and sleeve for accuracy in angulation
- Depth control with stop drills

FULLY GUIDED IMPLANT INSERTION

- Implant driver fits the sleeve, for a fully guided insertion with physical depth control
- Offset: 10 mm



1. DATA ACQUISITION
3D (CB)CT scan (DICOM)
Intraoral or lab scanning
(STL images)


3. SURGICAL GUIDE PRODUCTION
The surgical guide must contain
the sleeves that guide the
instruments and the implants.









2. VIRTUAL PLANNING
Implant positioned respecting the
patient's anatomy and prosthetic
outcome. Neodent® **EasyGuide** is
compatible with major software.

4. SURGICAL PROCEDURE
Neodent® **EasyGuide** presents
two surgical kits, selected
according to the implant
diameter.

EASYGUIDE KIT NARROW/REGULAR • Ø 3.5, Ø 3.75



EASYGUIDE KIT REGULAR/WIDE • Ø 4.0, Ø 4.3, Ø 5.0



Neodent® EasyGuide Kits

Neodent® EasyGuide Kit for Narrow/Regular Diameter Implants

Autoclavable polymer case.
The kit allows the installation of Helix GM® Implants of Ø3.5 and Ø3.75 in all bone types, using the Neodent® EasyGuide Guided Surgery Technique.



Articles

- | | | | | | |
|---------|--|---------|---------------------------------------|---------|----------------------------------|
| 110.313 | EasyGuide Kit Narrow/Reg. Diam. Tray | 103.551 | Narrow Tapered Drill D3.5/3.75X10 | 125.142 | Fixation Clamp - 3 units per kit |
| 125.170 | GM Narrow Stabilizer - 3 units per kit | 103.552 | Narrow Tapered Drill D3.5/3.75X11.5 | 129.034 | Depth Probe |
| 105.161 | GM Narrow Driver for Contra-angle | 103.553 | Narrow Tapered Drill D3.5/3.75X13 | 104.050 | Torque Wrench |
| 105.162 | GM Narrow Driver for Torque Wrench | 103.554 | Narrow Tapered Drill D3.75X8 | | |
| 103.583 | Narrow Mucosa Punch | 103.555 | Narrow Tapered Drill D3.75X10 | | |
| 103.519 | Narrow Bone Leveling Drill | 103.556 | Narrow Tapered Drill D3.75X11.5 | | |
| 103.545 | Narrow Initial Drill | 103.557 | Narrow Tapered Drill D3.75X13 | | |
| 103.546 | Narrow Tapered Drill D3.5X8 | 105.160 | Long Neo Screwdriver for Contra-angle | | |
| 103.547 | Narrow Tapered Drill D3.5X10 | 104.060 | Neo Manual Screwdriver (Medium) | | |
| 103.548 | Narrow Tapered Drill D3.5X11.5 | 103.558 | Drill for Palatal Setter | | |
| 103.549 | Narrow Tapered Drill D3.5X13 | 125.176 | Palatal Setter | | |
| 103.550 | Narrow Tapered Drill D3.5/3.75X8 | 103.395 | Guided Surgery Drill 1.3 | | |

Note: Items that compose Neodent® Kits are sold separately.

Neodent® EasyGuide Kit for Regular/Wide Diameter Implants

Autoclavable polymer case.
The kit allows the installation of Helix GM® Implants of Ø4.0, Ø4.3 and Ø5.0 in all bone types, using the Neodent® EasyGuide Guided Surgery Technique.



Articles

- | | | | | | |
|---------|---|---------|-------------------------------------|---------|---------------------------------------|
| 110.314 | EasyGuide Kit Reg./Wide Diam. Tray | 103.530 | Regular Tapered Drill D4.0X10 | 103.542 | Regular Tapered Drill D5.0X10 |
| 125.171 | GM Regular Stabilizer - 3 units per kit | 103.531 | Regular Tapered Drill D4.0X11.5 | 103.543 | Regular Tapered Drill D5.0X11.5 |
| 105.163 | GM Regular Driver for Contra-angle | 103.532 | Regular Tapered Drill D4.0X13 | 103.544 | Regular Tapered Drill D5.0X13 |
| 105.164 | GM Regular Driver for Torque Wrench | 103.533 | Regular Tapered Drill D4.0/4.3X8 | 105.160 | Long Neo Screwdriver for Contra-angle |
| 103.584 | Regular Mucosa Punch | 103.534 | Regular Tapered Drill D4.0/4.3X10 | 104.060 | Neo Manual Screwdriver (Medium) |
| 103.518 | Regular Bone Leveling Drill | 103.535 | Regular Tapered Drill D4.0/4.3X11.5 | 103.558 | Drill for Palatal Setter |
| 103.520 | Regular Initial Drill | 103.536 | Regular Tapered Drill D4.0/4.3X13 | 125.176 | Palatal Setter |
| 103.521 | Regular Tapered Drill D2.7X8 | 103.537 | Regular Tapered Drill D4.3/5.0X8 | 103.395 | Guided Surgery Drill 1.3 |
| 103.522 | Regular Tapered Drill D2.7X10 | 103.538 | Regular Tapered Drill D4.3/5.0X10 | 125.142 | Fixation Clamp - 3 units per kit |
| 103.523 | Regular Tapered Drill D2.7X11.5 | 103.539 | Regular Tapered Drill D4.3/5.0X11.5 | 129.034 | Depth Probe |
| 103.524 | Regular Tapered Drill D2.7X13 | 103.540 | Regular Tapered Drill D4.3/5.0X13 | 104.050 | Torque Wrench |
| 103.529 | Regular Tapered Drill D4.0X8 | 103.541 | Regular Tapered Drill D5.0X8 | | |

Note: Items that compose Neodent® Kits are sold separately.



Neodent® EasyGuide Instruments



Narrow Tapered Drills

- :: Available in surgical steel;
- :: For Helix GM® implants with Ø3.5 and Ø3.75 in diameter;
- :: Built-in titanium stops for a fully-guided procedure, matching the color of the sleeve of the surgical guide;
- :: Color code according to implant diameter;
- :: Laser-marked length.

	Ø 3.5	Ø 3.5/3.75	Ø 3.75
8.0	103.546	103.550	103.554
10.0	103.547	103.551	103.555
11.5	103.548	103.552	103.556
13.0	103.549	103.553	103.557



Regular Tapered Drills

- :: Available in surgical steel;
- :: For Helix GM® implants with Ø4.0, Ø4.3 and Ø5.0 in diameter;
- :: Built-in titanium stops for a fully-guided procedure matching the color of the sleeve of the surgical guide;
- :: Color code according to implant diameter;
- :: Laser-marked length.

	Ø 2.7	Ø 4.0	Ø 4.0/4.3	Ø 4.3/5.0	Ø 5.0
8.0	103.521	103.529	103.533	103.537	103.541
10.0	103.522	103.530	103.534	103.538	103.542
11.5	103.523	103.531	103.535	103.539	103.543
13.0	103.524	103.532	103.536	103.540	103.544



Guided Surgery Drill 1.3 and Guide Clamp

- :: Drill available in stainless steel;
- :: Guide Clamp available in titanium;
- :: For initial fixation of the surgical guide.

Drill Ø 1.3	Guide Clamp
103.395	125.142



Drill and Palatal Setter

- :: Drill and Palatal Setter available in stainless steel;
- :: Palatal Setter placed with the GM Implant Driver for Contra-angle;
- :: Maximum torque of 20 N.cm.

Drill	Palatal Setter
103.558	125.176



Mucosa Punches

- :: Available in stainless steel;
- :: To remove the mucosa before beginning the osteotomy.
- :: Rotation recommended: 60 rpm.

Narrow	Regular
103.583	103.584



Bone Leveling Drills

- :: Available in stainless steel;
- :: Built-in titanium stops matching the color of the sleeve of the surgical guide;
- :: For flattening bone surface before osteotomy.

Narrow	Regular
103.519	103.518



Initial Drills

- :: Available in stainless steel;
- :: Built-in titanium stops matching the color of the sleeve of the surgical guide;;
- :: For rupture of the cortical bone.

Narrow	Regular
103.545	103.520



GM Drivers for Contra-Angle

- :: Available in stainless steel;
- :: Color-coded according to the sleeve of the surgical guide;
- :: To start the implant placement through the surgical guide;
- :: Maximum torque 35 N.cm.

Narrow	Regular
105.161	105.163



GM Drivers for Torque Wrench

- :: Available in stainless steel;
- :: To finish the implant placement through the surgical guide;
- :: Maximum torque 60 N.cm.

Narrow	Regular
105.162	105.164



Guide Stabilizers

- :: Available in titanium;
- :: Color-coded according to the sleeve of the surgical guide;
- :: Additional fixation of the surgical guide.

Narrow	Regular
125.170	125.171

Depth Probe

- :: Available in titanium;
- :: With marks matching the Helix GM® implant lengths.



129.034



Neo Manual Screwdriver

- :: Available in surgical steel and titanium.

Medium
25 mm

104.060



Neo Screwdriver Torque Connection
- Contra-angle

- :: Available in stainless steel;
- :: Maximum torque 20 N.cm.

105.160



Torque Wrench

- :: Available in surgical steel;
- :: Fitting for square connections;
- :: Collapsible Wrench that allows for proper assembly and cleaning.

104.050

Sleeves for Neodent® EasyGuide

- :: Available in titanium;
- :: Sold in bags with 10 units each.



125.165 Regular Sleeve D5.2



125.168 Narrow Sleeve D3.93



125.177 Sleeve for Palatal Setter



125.143 Sleeve for Fixation Clamp

A SMILE FOR EVERYONE

NEODENT® NEOARCH® IMMEDIATE FIXED FULL-ARCH SOLUTION

Increasing expectations for shortened treatment duration represent a significant challenge for dental professionals especially in patients with anatomical deficiencies. The Neodent® Implant System offers an optimized solution for immediate fixed treatment protocols in edentulous patients even with severe atrophic maxilla. Neodent® NeoArch® allows to significantly improve patient satisfaction and quality of life by immediately restoring function and esthetics ⁽¹⁰⁾.



Immediate function resulting in shorter treatment times.

- Different implants techniques to avoid the use of grafting procedure⁽¹¹⁾.
- Optimized implant design to achieve high primary stability in all bone types⁽¹²⁾.



Immediate natural-looking esthetics with versatile restorative options.

- A broad gingival height abutment range to cater the patient's needs.
- Options of straight and angled abutments (17°, 30° and 45°).



Immediate peace of mind thanks to a stable foundation.

- One connection regardless of the diameters.
- Unique connection combining Platform Switching associated with a deep 16° Morse taper including an internal indexation.

SOLUTIONS FOR ALL CLINICAL NEEDS

A implant system designed for predictable immediate treatments in all bone types even with different conditions of the residual alveolar bone.



Helix GM®



Helix GM® Long



BONE RESORPTION

Helix GM® Long

PRODUCT FEATURES:

- Implants Description:
- Full dual tapered implant;
 - Hybrid contour with a cylindrical coronal part and conical on the apical area;
 - Active apex including a soft rounded small tip and helicoidal flutes;
 - Dynamic progressive thread design: from compressing trapezoidal threads on the coronal area to self-tapping threads on the apical part;
 - Double lead threaded implant;
 - Holder integrated to the implant body, which adapt in the packaging;
 - Neoporos surface;
 - Grand Morse® connection.

- Indications:
- Indicated for surgical intraoral installation, in bone types III/IV for cases of total or partial edentulism and for multiple-unit prostheses.


- Drilling features:
- For infraosseous positioning it is recommended to add 1 to 2 mm in length to the implant during surgical instrumentation.
 - Drilling speed: 500-800 rpm;
 - Implant insertion speed: 30 rpm;
 - Maximum torque for implant placement: 60 N.cm.

Available with:


NeoPoros®



Drill Sequence










	Initial 103.453	Ø 2.35 103.462	Ø 3.75 103.463	Ø 4.0 103.464
Ø 3.75 mm	Optional	✓	✓	
Ø 4.0 mm	Optional	✓	✓	✓

Bone types III and IV 

The procedure can be with Guided Surgery. Check the instruments for more information.

Helix GM® Long implants


	20.0 mm	22.5 mm	25.0 mm
Ø 3.75			
NeoPoros	109.1043	109.1044	109.1045
Ø 4.0			
NeoPoros	109.1046	109.1047	109.1048



GM Healing Abutment


Profile	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm
Ø 3.3	106.207	106.208	106.209	106.210	106.211	106.212
Ø 4.5	106.213	106.214	106.215	106.216	106.217	106.218

:: Use the manual Neo Screwdriver (104.060);
:: Do not exceed the insertion torque of 10 N.cm.



GM Customizable Healing Abutments

Profile	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm	6.5 mm
Ø 5.5	106.223	106.224	106.225	106.226	106.227	
Ø 7.0		106.228	106.229	106.230	106.231	106.232





GM Cover Screw


	0 mm	2 mm
	117.021	117.022

:: Use the manual Neo Screwdriver (104.060);
:: Do not exceed the insertion torque of 10 N.cm.

GM Mini Conical Abutment


 Multiple-unit screw-retained prosthesis

 Ø 4.8 mm



Measurements GM Mini Conical Abutment

Installation Sequence




Drivers


 Hexagonal Prosthetic Driver +  Torque Wrench

 Neo Screwdriver Torque Connection +  Torque Wrench

 Neo Screwdriver Torque Connection +  Manual Screwdriver Torque

Accessories

 Mini Conical Abutment Polishing Protector 123.008

 Replacement Coping Screw
116.269 Titanium
116.270 Neotorque*

*Application of a film carbon-based coat that provides a lower friction coefficient, resulting in increased pre-load.

*The 45° Mini Conical Abutment is indicated for use only with Helix GM® Long and Zygoma GM™.

NeoArch® Kits

Helix GM® Long Compact Surgical Kit

Autoclavable polymer case.



Articles

110.300	Helix GM® Long Compact Surgical Kit Case	103.453	Helix GM® Long Initial Drill 2.0mm	105.143	Regular Guided Surgery GM Connection for Torque Wrench
103.395	Guided Surgery Drill 1.3mm	103.462	Twist Drill For Helix GM® Long 2.35mm	105.140	Regular Guided Surgery GM Connection - Contra-angle
125.100	Guided Surgery Guide Clamp	103.463	Twist Drill For Helix GM® Long 3.75mm	104.060	Neo Manual Screwdriver (medium)
125.140	Drill Guide For NGS Helix GM® Long 2.0/2.35mm	103.464	Twist Drill For Helix GM® Long 4.0mm	105.129	GM Implant Driver - Torque Wrench (short)
125.141	Drill Guide For NGS Helix GM® Long 3.75/4.0mm	129.021	Helix GM® Long X-ray Positioner	105.131	GM Implant Driver - Contra-angle
103.459	Twist Drill For NGS Helix GM® Long 2.35mm	128.032	GM Angle Measurer 17°	104.050	Torque Wrench
103.460	Twist Drill For NGS Helix GM® Long 3.75mm	128.033	GM Angle Measurer 30°		
103.461	Twist Drill For NGS Helix GM® Long 4.0mm	128.034	GM Angle Measurer 45°		

Note: Items that compose Neodent® Kits are sold separately.

NeoArch® Instruments



Helix GM® Long Drills

- :: Available in surgical steel;
- :: Drill sequence for Helix GM® Long implants.


Initial	Ø 2.35	Ø 3.75	Ø 4.0
103.453	103.462	103.463	103.464



Helix GM® Long Drills for Guided Surgery

- :: Available in surgical steel;
- :: Drill sequence for Helix GM® Long implants on Guided Surgery.

Ø 2.35	Ø 3.75	Ø 4.0
103.459	103.460	103.461



GM Height Measurer

- :: Available in titanium;
- :: For selecting GM prosthetic abutments;
- :: Marks corresponding to transmucosa heights.
- :: Can be used as X-Ray Positioner.

128.028



GM Implant Driver - Contra-Angle

- :: To capture the implant directly from the packaging;
- :: To place GM Implants with contra-angle, or attached to a manual driver for contra-angle connections (104.028) for hand placement;
- :: With six dimples to indicate the hex index face position;
- :: The laser marks indicate the depth of implant placement, bone level, 1 and 2mm infra-bone and last marking (3mm) biological space;
- :: Maximum torque 35 N.cm.

105.131



GM Implant Driver - Torque Wrench

- :: To place GM Implants with the Torque Wrench (104.050);
- :: With six marks to indicate the hex index face position;
- :: The laser marks indicate the depth of implant placement, bone level, 1 and 2mm infra-bone and last marking (3mm) biological space;
- :: Maximum torque: 60 N.cm.

Short 22 mm	Long 30 mm
105.129	105.130



Neo Screwdriver Torque Connection - Torque Wrench

- :: Available in surgical steel;
- :: Yellow color for line identification.

Short 16.5 mm	Medium 22 mm	Long 32 mm
105.133	105.132	105.157



Neo Manual Screwdriver

- :: Available in surgical steel;
- :: Yellow color for line identification.

Short 21 mm	Medium 25 mm	Long 37 mm
104.058	104.060	104.072



Neo Screwdriver Torque Connection - Contra-angle

- :: Available in surgical steel;
- :: Yellow color for line identification;
- :: Medium Neo Screwdriver Torque Connection
- :: Extra Short Neo Screwdriver Torque Connection - Contra-angle (105.146) recommended for Impression Copings, Cover Screws and Healing Abutments.

Extra Short 16.5 mm	Short 24 mm	Long 31 mm
105.146	105.135	105.160



Hexagonal Prosthetic Driver

- :: Available in surgical steel;
- :: To install and apply torque over straight GM Mini Conical Abutments and GM Micro Abutments;
- :: Yellow color for line identification.

- Contra-angle
- 105.138
- Torque Wrench
- 105.137



GM Bone Profile Drill with Guide

- :: Available in surgical steel;
- :: Used in the surgical second step;
- :: Conforms the bone around the implant platform, preparing the emergence profile to be suitable to prosthetic components.

103.424



GM Angle Measurer

- :: Available in titanium;
- :: Angles: 17°, 30° and 45°;
- :: To a more accurate selection and planning of the abutments angulation during the prosthetic phase.

- 17°
- 128.032
- 30°
- 128.033
- 45°
- 128.034



Helix GM® Long Drill Guide for Guided Surgery

- :: Instrument with the purpose of guiding the drills during the bone bed preparation according to the guided surgery technique.

- Ø 2.0/2.35
- 125.140
- Ø 3.75/4.0
- 125.141



Guided Surgery Drill 1.3 and Guide Clamp

- :: Drill available in surgical steel;
- :: Guide Clamp available in titanium;
- :: For initial fixation of the surgical guide.

- Drill Ø 1.3
- 103.395
- Guide Clamp
- 125.100



Guided Surgery GM Connection - Contra-Angle

- :: Available in stainless steel;
- :: To start the implant placement through the surgical guide.

Regular
105.140



Guided Surgery GM Connection - Torque Wrench

- :: Available in stainless steel;
- :: To finish the implant placement through the surgical guide.

Regular
105.143



Helix GM® Long X-ray Positioner

- :: Indicated for evaluation of the osteotomy depth in the implant placement procedure.

129.021



Torque Wrench

- :: Available in surgical steel;
- :: Fitting for square connections;
- :: Collapsible Wrench that allows for proper assembly cleaning;
- :: For full instructions see page 80.

104.050



Remover for Abutments with internal threads

- :: Available in surgical steel;
- :: To remove abutments with internal threads from the implants, after removal of the screws;
- :: Compatible with abutments with Neo removable Screws

- Long
- 130.118
- 130.114



Remover for Neo Screws

- :: Available in surgical steel;
- :: Compatible with Neo removable screws for abutments

- Long
- 130.119
- 130.115



- Long
- 130.117
- 130.116

Removal Sets for Abutments with internal threads and Neo Screws

- :: Available in surgical steel;
- :: To remove Neo Removable Screws and abutments with internal threads from the implants, after removal of the screws;
- :: Compatible with abutments with Neo removable Screws

GRAND MORSE® NEODENT® GUIDED SURGERY. GRAND POSSIBILITIES WITH A LIMITLESS SOLUTION

Patients’ expectations regarding tooth replacement are increasing and are even higher when it comes to treatment duration and esthetic outcomes. The Neodent® Guided Surgery helps clinicians to provide prosthetically driven treatments, enabling them to perform immediate protocols with peace of mind, fulfilling patients’ expectations.



DIFFERENTIATE YOUR PRACTICE WITH GUIDED SURGERY.



- Improve patient quality of life.**
- Functional with an immediate fixed restoration.
 - Esthetical with a personalized restoration and less bone remodeling ⁽¹³⁾.
 - Comfort by the reduction of operative and postoperative discomfort (e.g. reduced patient chair time).



- Access to more treatment options.**
- Reliable access to flapless surgery ⁽¹⁴⁻¹⁶⁾.
 - Designed to reduce bone grafting procedures.
 - Predictable immediate protocols.



- Increase patient acceptance.**
- Better communication building trust with patients.
 - Reliable treatment estimates from root to tooth including components and procedures.

SURGICAL PREDICTABILITY AND EFFICIENCY WITH A LIMITLESS SOLUTION.

Guided surgery is designed to reduce chair time and postoperative discomfort. It helps increasing implant positioning accuracy ⁽¹⁷⁾.




Complete
Helix® and Drive GM®
Implants portfolio



Convenient
Color-coded instruments
and symbol-marked



Flexible
2 sleeve height positions



Neodent® Guided Surgery Kit

Grand Morse® Guided Surgery Surgical Kit

Autoclavable polymer case.
The Kit allows the use of Helix GM® and Drive GM® Implants in the Guided Surgery technique.



Articles

- 110.296 GM Guided Surgery Surgical Kit Case
- 103.395 Guided Surgery 1.3
- 125.100 Guided Surgery Guide Clamp
- 103.429 Narrow Guided Surgery Punch - Contra-Angle
- 103.430 Regular Guided Surgery Punch - Contra-Angle
- 103.431 Wide Guided Surgery Punch - Contra-Angle
- 103.432 Guided Surgery Drill 2.0
- 103.433 Tapered Guided Surgery Drill 3.5*
- 103.434 Tapered Guided Surgery Drill 3.75*
- 103.435 Tapered Guided Surgery Drill 4.0*
- 103.436 Tapered Guided Surgery Drill 4.3*
- 103.437 Tapered Guided Surgery Drill 5.0*
- 103.438 Tapered Guided Surgery Drill 6.0*
- 105.139 Narrow Guided Surgery GM Connection - Contra-angle
- 105.140 Regular Guided Surgery GM Connection - Contra-angle
- 105.141 Wide Guided Surgery GM Connection - Contra-angle
- 105.142 Narrow Guided Surgery GM Connection for Torque Wrench
- 105.143 Regular Guided Surgery GM Connection for Torque Wrench
- 105.144 Wide Guided Surgery GM Connection for Torque Wrench
- 125.130 Narrow Guided Surgery GM Guide Stabilizer
- 125.131 Regular Guided Surgery GM Guide Stabilizer
- 125.132 Wide Guided Surgery GM Guide Stabilizer
- 125.133 Narrow Guided Surgery GM Guide Stabilizer (Long)
- 125.134 Regular Guided Surgery GM Guide Stabilizer (Long)
- 105.145 Guided Surgery GM H11 Connection for Torque Wrench
- 105.136 Neo Screwdriver Torque Connection - Contra-angle (Medium)

- 104.060 Neo Manual Screwdriver (Medium)
- 103.439 Tapered Contour Guided Surgery Drill 3.5*
- 103.440 Tapered Contour Guided Surgery Drill 3.75*
- 103.441 Tapered Contour Guided Surgery Drill 4.0*
- 103.442 Tapered Contour Guided Surgery Drill 4.3*
- 103.443 Tapered Contour Guided Surgery Drill 5.0*
- 103.444 Narrow Guided Surgery GM Pilot Drill 3.5
- 103.445 Regular Guided Surgery GM Pilot Drill 3.5
- 103.446 Guided Surgery GM Pilot Drill 3.75
- 103.447 Guided Surgery GM Pilot Drill 4.0
- 103.448 Guided Surgery GM Pilot Drill 4.3
- 103.449 Guided Surgery GM Pilot Drill 5.0
- 125.119 Narrow Guided Surgery Drill Guide 2.0/3.5
- 125.121 Regular Guided Surgery Drill Guide 2.0/3.5
- 125.122 Regular Guided Surgery Drill Guide 3.75/4.0
- 125.123 Regular Guided Surgery Drill Guide 4.3
- 125.126 Wide Guided Surgery Drill Guide 2.0/3.5
- 125.127 Wide Guided Surgery Drill Guide 4.0/4.3
- 125.128 Wide Guided Surgery Drill Guide 5.0/6.0
- 125.120 Narrow Tapered Contour Guided Surgery Drill Guide 3.5
- 125.124 Regular Tapered Contour Guided Surgery Drill Guide 3.5/3.75
- 125.125 Regular Tapered Contour Guided Surgery Drill Guide 4.0/4.3
- 125.129 Wide Tapered Contour Guided Surgery Drill Guide 5.0
- 129.001 Titanium Tweezers
- 104.050 Torque Wrench

Note: Items that compose Neodent® Kits are sold separately.
*Conventional guided surgery drills that can be replaced by the respective short version.

Neodent®

Guided Surgery Instruments



Guided Surgery Tapered Drills

:: Available in surgical steel;
:: Drill sequence for Helix GM® and Drive GM®
Implants in the guided surgery technique;
:: Fully guided technique with Short Drills indicated
for 8, 10 or 11.5 mm long implants.

Short 36.5 mm	Ø 2.0	Ø 3.5	Ø 3.75	Ø 4.0	Ø 4.3	Ø 5.0	Ø 6.0
	103.475	103.476	103.477	103.478	103.479	103.480	103.481
Regular 41 mm	103.432	103.433	103.434	103.435	103.436	103.437	103.438



Guided Surgery Drill 1.3 and Guide Clamp

:: Drill available in surgical steel;
:: Guide Clamp available in titanium;
:: For initial fixation of the surgical guide.

Drill Ø 1.3	Guide Clamp
103.395	125.100



Guided Surgery Tapered Contour Drills

:: Available in surgical steel;
:: Drill sequence for Helix GM® Implants in the guided
surgery technique for bone types I or II;
:: Fully guided technique with Short Drills indicated
for 8, 10 or 11.5 mm long implants.

Short 36.5 mm	Ø 3.5+	Ø 3.75+	Ø 4.0+	Ø 4.3+	Ø 5.0+
	103.482	103.483	103.484	103.485	103.486
Regular 41 mm	103.439	103.440	103.441	103.442	103.443



Guided Surgery Punch - Contra-Angle

:: Available in titanium;
:: Color-coded according to the sleeve
diameter;
:: To remove the mucosa before beginning
the osteotomy.

Narrow	Regular	Wide
103.429	103.430	103.431



Guided Surgery GM Pilot Drills

:: Available in surgical steel;
:: Color-coded according to the sleeve diameter;
:: Recommended for Helix GM® in bone types I or II;
:: Optional Drive GM® in bone types III or IV.

Narrow	Ø 3.5	103.444	Regular	Ø 3.5	103.445	Wide	Ø 5.0	103.449



Guided Surgery Drill Guides

:: Available in titanium and stainless steel;
:: Color-coded according to the sleeve diameter;
:: To fit in the sleeve in the surgical guide;
:: To be used with correspondent drill diameter
and type.

Narrow	Ø 2.0/3.5	125.119	Regular	Ø 2.0/3.5	125.121	Wide	Ø 2.0/3.5	125.126
Ø 3.5+	125.120	Ø 3.75/4.0	125.122	Ø 4.0/4.3	125.127	Ø 5.0/6.0	125.128	125.129
Ø 3.5+/3.75+	125.124	Ø 4.0+/4.3+	125.125	Ø 4.3	125.123	Ø 5.0+/5.5+	125.124	125.125



Guided Surgery GM Connection
- Contra-Angle

- :: Available in stainless steel;
- :: Color-coded according to the sleeve diameter;
- :: To start the implant placement through the surgical guide.

Narrow	Regular	Wide
105.139	105.140	105.141



Guided Surgery GM Connection
- Torque Wrench

- :: Available in stainless steel;
- :: Color-coded according to the sleeve diameter;
- :: To finish the implant placement through the surgical guide.

Narrow	Regular	Wide
105.142	105.143	105.144



Guided Surgery GM H 11 Connection
- Torque Wrench

- :: Available in stainless steel;
- :: To finish the implant placement through the surgical guide;
- :: To be used when the H11 sleeve height is chosen.

105.145



Guided Surgery Guide
Stabilizers

- :: Available in titanium;
- :: Color-coded according to the sleeve diameter;
- :: Additional fixation of the surgical guide.

Narrow	Regular	Wide
125.130	125.131	125.132



Guided Surgery Guide Stabilizers - Long

- :: Available in titanium;
- :: Additional fixation of the surgical guide;
- :: To be used when the H11 sleeve height is chosen.

Narrow	Regular
125.133	125.134

Sleeves for Neodent® Guided Surgery System


- :: Available in titanium;
- :: Sold in bags with 10 units each.



125.135	Sleeve for Narrow Guided Surgery System
125.136	Sleeve for Regular Guided Surgery System
125.137	Sleeve for Wide Guided Surgery System
125.138	Sleeve of Setter for Guided Surgery System

One Step Hybrid Technique

Technique that allows passive fitting, with no need for welding as the titanium coping is cemented to the substructure.
Used for multiple prostheses and reduces laboratory work times.



Neo Mini Conical Abutment One Step Hybrid Copings

:: For installation, use the Neo Torque Connection (105.132);
:: For torque control, use Torque Wrench (104.050).

Burn-out	Brass	Titanium
118.340	118.331	118.330



Neo Micro Conical Abutment One Step Hybrid Copings

:: For installation, use the Neo Torque Connection (105.132);
:: For torque control, use Torque Wrench (104.050).

Burn-out	Brass	Titanium
118.341	118.333	118.332



Neo Working Screw One Step Hybrid

:: For laboratory use.

116.271


Demonstration Sequence



Regularize the alveolar ridge.



Surgical drilling completed, obtaining adequate distance from distal implant in relation to the mental foramen with 7 mm Space Planning Instrument.



Placement of 4 Neodent® implants, according to their indication.



Placement of corresponding Neodent® Abutments.



Placement of Impression Copings, splinted with acrylic resin.



Positioning of Multifunctional Guide to obtain intermaxillary correlation. Soft silicone is injected to take the soft tissue impression.



Removal of Multi-Functional Guide and placement of Analogs to the impression copings.



Working model with artificial gum.



Burn-out One Step Hybrid Coping, Brass One Step Hybrid Coping, grooved Titanium One Step Hybrid Coping. The last one with lower dimensions than the brass one, which compensates using the mill.



Brass Copings are placed over analogs, then Burn-out Copings are fixed by working screws.



Castable ring with waxed framework.



Cast framework.



Place the framework over the stone model.



Please note cementing area.



Cementing with Panavia the structure over the titanium copings.



Final inside-mouth view.

Distal Bar Technique

Technique used to ease mandible rehabilitation, through a provisional hybrid type prostheses supported by implants.



Neo Distal Bar Coping

- :: Available in titanium;
- :: Retainers to ease joining with acrylic resin;
- :: Recommended torque: 10 N.cm;
- :: For torque, use Neo Screwdriver (105.132)

118.308

Neo Distal Bar

- :: Recommended for distal Implants to reinforce the cantilever.

125.116

Polishing Protector

- :: Available in surgical steel;
- :: Protection for the lab polishing.

123.008

Demonstration Sequence



1 Neodent® Abutments placed.



2 Prosthesis wearing, keeping posterior region integrity.




3 Place the copings into the central Implants and Distal Bar to distal Implants.



4 Proof of inferior prostheses wearing (centered occlusion position, no interference on copings).



5 Placement of rubber dam over copings to protect soft tissues.



6 Apply selfpolymerizing acrylic resin on and between the copings.



7 Apply to worn area in lower prosthesis, repositioning inside mouth. Keep patient in occlusion until total polymerization.



8 Remove the inferior prosthesis after resin is polymerized. Copings already captured.



9 Adjustments, finishing and polishing procedures of inferior prosthesis with polishing protectors.



10 Placed provisional implant supported prosthesis.



11 Final inside-mouth posterior view.



Visit www.neodent.com/cadcam to download the digital files to work with Neodent® Titanium Bases, Titanium Blocks, Abutments, Mini Conical Abutments, Micro Abutments, Universal Abutments, One Step Hybrid Copings, Scanbodies and Hybrid Repositionable Analogs. Libraries are available for the following companies: exocad GmbH, Amann Girrbach AG Inc, Dental Wings Inc and 3Shape A/S.

Scanbody

Neodent® Scanbodies can be used for scanning and digitalization of the patient or model providing accuracy in determining the analog position.



- 108.183 GM Exact Implant Intraoral Scanbody
- 108.181 GM Exact Implant Scanbody (for model)
- 108.196 GM Mini Conical Abutment Scanbody (intraoral and model)
- 108.197 GM Micro Abutment (intraoral and model)
- 108.198 GM Abutment (intraoral and model)



Compatible with Neo Screwdriver

Hybrid Repositionable Analog

Neodent® Hybrid Repositionable Analogs can be used in prototyped models, produced by 3D printers, or conventional plaster models.



- 101.103 GM Hybrid Repositionable Analog 3.5/3.75
- 101.089 GM Hybrid Repositionable Analog 4.0/4.3
- 101.090 GM Hybrid Repositionable Analog 5.0/6.0
- 101.091 Micro Abutment Hybrid Repositionable Analog
- 101.092 Mini Conical Abutment Hybrid Repositionable Analog
- 101.097 Universal Abutment Hybrid Repositionable Analog 3.3X4
- 101.098 Universal Abutment Hybrid Repositionable Analog 3.3X6
- 101.099 Universal Abutment Hybrid Repositionable Analog 4.5X4
- 101.100 Universal Abutment Hybrid Repositionable Analog 4.5X6
- 101.101 GM Abutment Hybrid Repositionable Analog

Torque Wrench

- :: Available in surgical steel;
- :: Extremely safe (lower than 5% variation);
- :: Fitting for square connections;
- :: Collapsible Wrench that allows for proper assembly cleaning.


104.050



Operational Instructions

The Neodent® Torque Wrench was designed to allow the necessary torque to be applied and simultaneous verification of that torque with the same Instrument.

All that is needed is to apply force to the wrench handle **1** (never the wrench body) until the value marked on the LATERAL SCALE **2** corresponds to the desired torque.



The wrench function works in both directions, by simply pulling and turning the driver's pin 180°. However, the torque measurements work only lockwise.

•WARNING: When inverting the torque direction, the gear may come loose from the driver body and fall. Therefore, this inversion should only be done with the driver connected to a part or outside the patient's mouth.



The Neodent® Torque Wrench comes with pre-calibrated torques



Titanium Tweezers

- :: To handle implants;
- :: New Tweezer system that prevents deviation in the active bit;
- :: Millimeter scale for checking during procedures;
- :: Self-locking implant.

129.001



Depth Probe

- :: Available in titanium;
- :: To probe preparations and analyze depth;
- :: Millimeter scale for checking during procedures.

129.004



7 and 9 mm Space Planning Instrument

- :: Available in surgical steel;
- :: Recommended for prosthetic/surgical planning.
- :: 7 and 9 mm marks.

128.026



Surgical Labial Retractor

- :: Available in surgical steel;
- :: Rounded edges to minimize surgical trauma.

124.001



Columbia Retractor

- :: Available in surgical steel;
- :: Rounded edges to minimize surgical trauma.

124.003



Scapel Handle

- :: Available in surgical steel;
- :: For standard scalpel blade use;
- :: Blade not included.

129.008




Bivers Handle

- :: Available in surgical steel;
- :: Non-traumatic extraction for implant placement;
- :: Similar to a periostome.

129.002





17 mm

13 mm

9 mm

15 mm

11 mm

7 mm

1.8 mm

2.5 mm

3.0 mm

3.5 mm

4.0 mm

4.5 mm

110.154

110.155

110.156

110.157

110.158

110.159

Concave Osteotome


:: Available in surgical steel;

:: Concave active cutting bit for nontraumatic lifting the floor of the maxillary sinus;

:: Used to prepare the surgical alveolus for Implant placement in the posterior maxillary region with low bone height;

:: Marks from 7 to 17mm.

:: Marks from 7 to 17mm.



17 mm

13 mm

9 mm

15 mm

11 mm

7 mm

1.8 mm

2.5 mm

3.0 mm

3.5 mm

110.160

110.161

110.162

110.163

Convex Osteotome

:: Available in surgical steel;

:: Convex active bit;

:: Used when the bone width is insufficient, demanding bone compression and expansion before placing the implant;

:: Marks from 7 to 17mm.



110.262

Osteotomes Kit Case

:: Available in polymer;

:: Autoclavable;

:: Osteotomes sold separately.



126.001



Surgical Hammer

:: Available in surgical steel;

:: Polymer active bit;

:: Used in compactors and expanders;

:: Weight: 130g.



0.35 mm

Ø 3.3

Ø 3.5

Ø 3.75

Ø 4.1

Ø 4.3

Ø 5.0

Ø 8.0

103.051

103.490

103.491

103.026

103.087

103.027

103.028

Trephine Bur

:: Available in surgical steel;

:: Collecting bone cylinder;

:: Implant removal.



1

3

4

5

7

126.008

126.009

126.010

126.011

126.012

Sinus Lift Curette

:: Available in surgical steel;

:: Used to displace the Sinusal Membrane.




110.270

Complement Case

:: Available in autoclavable polymer;

:: Used to organize drills and auxilliary connections.




104.047

Handle Implant Driver

:: Available in stainless steel;


:: Manual implant placement.



104.036

Analog Handle

:: Used for tightening analogs and milling prosthetic abutments.



Guide

Pin

103.092

103.093

Prosthetic Surgical Guide

:: Available in titanium;

:: Abutments to prepare the surgical guide;

:: Prosthetic guide inner diameter 2 mm

:: Heights 6 and 10 mm;

:: Surgical Guide: package with 10 units (5 units of 10 mm and 5 units of 6 mm);

:: Surgical Guide Pin: package with 5 units

Najlepsze decyzje są oparte na faktach.

<div>25 LATA</div> <div>Dostarczamy wysokiej jakości rozwiązania implantologiczne od ponad 25 lat</div>	<div>NR2</div> <div>Jesteśmy drugą co do wielkości firmą implantologiczną na świecie</div>	<div>Więcej niż 270,000</div> <div>pacjentów było leczonych implantami Neodent w Brazylii w 2014 roku</div>
<div>Udokumentowany wskaźnik powodzenia</div> <div>99.7%</div> <div>u 3578 pacjentów</div>	<div>Więcej niż 250+</div> <div>naukowych artykułów i publikacji</div>	<div>NR1</div> <div>Neodent jest liderem w Ameryce Łacińskiej - drugim co do wielkości rynku zbytu implantów.</div>
<div>100%</div> <div>wsparcie dla naszych klientów</div>	<div>45,000</div> <div>profesjonalistów na świecie używa implantów Neodent</div>	<div>1,600,000</div> <div>Sprzedajemy corocznie ponad milion implantów</div>
<div>Edukacja, nauka, wsparcie i badania naukowe w Curitiba, Brazylia</div> <div>ILAPEO</div>	<div>Straumann zwiększył swoje udziały w Neodent do 100%</div> <div>(w kwietniu 2015 roku)</div>	<div>Ponad 3,000</div> <div>lekarzy zostało przeszkolonych w ośrodkach naukowych Neodent w Brazylii</div>

ILAPEO

Neodent kładzie dużą wagę na edukację i szkolenia w celu zapewnienia jakości oraz standardów leczenia i opieki nad pacjentem. Czyni to poprzez ILAPEO (Latin American Institute of Dental Education), jeden z największych na świecie ośrodków szkolenia dedykowany implantologii, który oferuje kursy dla tysięcy uczestników z całego świata każdego roku. ILAPEO proponuje kursy podyplomowe i kursy mistrzowskie dla klinicystów. Jest to wspaniały powód, aby odwiedzić Brazylię....

NEODENT - DOŻYWOTNIA GWARANCJA

Celem programu Neodent Originals jest zapewnienie jakości i trwałości produktów Neodent. Dlatego stomatolodzy, którzy używają w swoich procedurach produktów Neodent mogą liczyć na gwarancję dotyczącą każdego rozwiązania.

IMPLANTY: Dożywotnia gwarancja. Wymiana na taki sam lub podobny implant.

KOMPONENTY IMPLANTÓW: 10 lat gwarancji na łączniki metalowe.



Znakomity marketing Neodent usprawni Twoją komunikację z pacjentami.

Odwiedź nową stronę dla pacjentów Neodent:

www.mojimplant.com.pl