

Straumann® CARES® X-Stream™

More than efficient.
X-Streamly time-saving.




 **straumann**
simply doing more

Straumann® CARES® X-Stream™

WHAT IS CARES® X-STREAM™?

CARES® X-Stream™ is a highly efficient CAD/CAM workflow that allows you to gain time by creating both, the abutment and the corresponding restoration (crown, bridge or over-denture) in a single streamlined process. Design the restoration first and let the software design what's beneath it. Both parts are then milled simultaneously and shipped in one package.



YOUR BENEFITS

- **Efficiency:** You want to step up the productivity and efficiency of your business? Save up to 50% of the steps compared to conventional CAD/CAM processing.
- **Precision:** Fed up with checking the fit? Our validated milling center creates the right fit for you.
- **Flexibility:** Monolithic or framework? You choose the abutments, materials and techniques that match your professional needs.
- **Connectivity:** No Straumann scanner? Use one of our partner systems or benefit from the CARES® Scan & Shape¹ service.
- **Quality:** The Straumann® original connection ensures a perfect fit with Straumann implants.

¹ The CARES® Scan & Shape service may not be available in your country, please check with your local sales representative for further information.

Product overview

MATERIALS



3M™ ESPE™ Lava™ Plus
Outstanding esthetics with true colors and a plus in translucency



zerion®
Low translucent zirconium-dioxide for highly stable frameworks



IPS e.max® CAD
All ceramic restorations for efficient esthetics



coron®
Cobalt-chrome frameworks: Veneering process comparable to noble casting alloys



ticon®
Highly reliable and precise titanium frameworks



polycon® ae²
Full-contour temporary PMMA restorations for high efficiency

ABUTMENTS



CARES® Abutment, ceramic (ZrO₂)
For highly esthetic & customized single and multi-tooth restorations



Variobase® for crown
For efficient single-tooth restorations



CARES® Abutment, titanium
For customized single and multi-tooth restorations



Variobase® for bridge/bar
For efficient multi-tooth restorations



Variobase® coping for bridge/bar
For efficient multi-tooth restorations on screw-retained abutments

² Can be used in combination with all abutments (for single- and multi-tooth restorations)

International Headquarters

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

www.straumann.com

3M™, ESPE™ and Lava™ are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. IPS e.max® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.

© Institut Straumann AG, 2015. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.