




Straumann® Modular Cassette Selection Guide

Content

A Module Full List	5
Standard Plus, Standard Plus Short, Bone Level (SP, SPS, BL)	10
Tapered Effect (TE)	14
Straumann® PURE Ceramic Implant System (CIM/CI)	18
Bone Level Tapered (BLT)	24
Bone Level X (BLX)	28
Bone Level X (BLX) guided	32

Purpose of this Guide

The selection guide helps to define the required Modular Cassette parts and instruments needed for a specific indication. Each chapter includes an order list as well as a set-up overview which shows how to assemble the parts.



A Module Full List


This chapter shows all mandatory *and* optional instruments for the A Module.

A Module



A Module		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module

Ratchet Tray



Ratchet Tray			
Art. No.	Image	Product	Indication
041.766		A Modul Ratchet Tray	
046.119		Ratchet	
046.049		Torque Control Device	
066.1100		Torque Control Device	BLX
046.064		Holding Key	

Released Aid Tray




Released Aid Tray			
Art. No.	Image	Product	Indication
041.767		Released Aid Tray	
026.0073		Released Aid S for Loxim	BLT SC
026.2558		Released Aid N for Loxim	S, SP, TE, BL, BLT
026.4558		Released Aid R/W for Loxim	S, SP, TE, BL, BLT

Stop Key Tray



Stop Key Tray		
Art. No.	Image	Product
041.768		Stop Key Tray
034.006		Stop Key, guided

Tweezer Tray




Tweezer Tray		
Art. No.	Image	Product
041.770		Tweezer Tray
046.110		Instrument Tweezers
046.095		Tweezers with t.p.c. tips

Planning Aid Tray



Planning Aid Tray			
Art. No.	Image	Product	Indication
041.769		Planning Aid Tray	S, SP, TE, CIM/CI
046.147		Diagnostic T	S, SP, TE, CIM/CI, BL, BLT
046.148		Tissue Level Implant Distance Indicator	S, SP, TE, CIM/CI
026.0901		Bone Level Implant Distance Indicator	BL, BLT

Implant Depth Gauge Tray



Implant Depth Gauge Tray		
Art. No.	Image	Product
041.771		Implant Depth Gauge Tray
066.2000		Implant Depth Gauge

SCS Screwdriver Set



SCS Screwdriver Set		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS screwdriver for Ratchet, extra short
046.401		SCS screwdriver for Ratchet, short
046.402		SCS screwdriver for Ratchet, long
046.410		SCS screwdriver for Handpiece, extra short
046.411		SCS screwdriver for Handpiece, short
046.412		SCS screwdriver for Handpiece, long

AS Screwdriver Set



AS Screwdriver Set		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
046.789		AS Screwdriver for Handpiece, extra short
046.790		AS Screwdriver for Handpiece, short
046.791		AS Screwdriver for Handpiece, long
046.786		AS Screwdriver for Ratchet, extra short*
046.787		AS Screwdriver for Ratchet, short*
046.788		AS Screwdriver for Ratchet, long*
041.765		Grommet Tray 3 medium + 3 large

* Store the handling aid of the AS Screwdriver for Ratchet separately in the 041.765 Grommet tray

Bone Profiler Set



Bone Profiler Set			
Art. No.	Image	Product	Indication
041.762		Grommet Tray 6 small	
044.304		Bone Profiler RN	S, SP, TE
044.305		Bone Profiler WN	S, SP, TE
026.0070		Bone Profiler, oval Platform SC	BL, BLT
026.0022		BL Bone Profiler 1	BL, BLT, BLX
026.0023		BL Bone Profiler 2	BL, BLT, BLX
026.0024		BL Bone Profiler 3	BL, BLT, BLX

48h Tool



48h Tool			
Art. No.	Image	Product	Indication
041.763		Grommet Tray 6 large	
80103		48 hour Explantation Device RN/WN	S, SP, TE
80087		48 hour Explantation Device RN/WN	S, SP, TE
031.080		48 hour Explantation Device, Straumann® PURE Ceramic Implant (Monotype) RD	CIM/CI
031.081		48 hour Explantation Device, Straumann® PURE Ceramic Implant (Monotype) ND	CIM/CI
032.053		48h Explantation Device, Straumann® PURE Ceramic Implant RD	CIM/CI
026.0072		48 hour Explantation Device SC	BLT
026.2048		48 hour Explantation Device NC	BL, BLT
026.4048		48 hour Explantation Device RC	BL, BLT

* Note: CIM/CI PURE 48h explantation devices are for single use only.



Order List: Standard Plus, Standard Plus Short, Bone Level

A Module – Order list

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device
046.064		Holding Key
041.767		Released Aid Tray
026.2558		Released Aid N for Loxim
026.4558		Released Aid R/W for Loxim
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS screwdriver for Ratchet, extra-short
046.401		SCS screwdriver for Ratchet, short
046.402		SCS screwdriver for Ratchet, long
046.410		SCS screwdriver for Handpiece, extra short
046.411		SCS screwdriver for Handpiece, short
046.412		SCS screwdriver for Handpiece, long

A Module Tissue Level Specific Tools

Add for Standard Plus and Standard Plus Short Implant line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
044.304		Bone Profiler RN
044.305		Bone Profiler WN
80103		48 hour Explantation Device RN/WN
80087		48 hour Explantation Device RN/WN

A Module Bone Level Specific Tools

Add for Bone Level Implant Line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3
026.2048		48 hour Explantation Device NC
026.4048		48 hour Explantation Device RC

For additional tools (planning aid, tweezers, AS Screwdriver), please refer to page 5 *A Module full list*.



Order List: Standard Plus, Standard Plus Short, Bone Level

B Module – Order list

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module
041.779		B Module, Parallel Walled Tray
040.563		Drill extender
044.003		Round bur Ø 2.3 mm
044.004		Round bur Ø 3.1 mm
026.0056		Needle drill long Ø 1.6 mm
044.210		Pilot Drill 1, short, Ø 2.2 mm
044.214		Pilot Drill 2, short, Ø 2.8 mm
044.250		Twist Drill PRO, short, Ø 3.5 mm
044.254		Twist Drill PRO, short, Ø 4.2 mm
044.211		Pilot Drill 1, long, Ø 2.2 mm
044.215		Pilot Drill 2, long, Ø 2.8 mm
044.251		Twist Drill PRO, long, Ø 3.5 mm
044.255		Twist Drill PRO, long, Ø 4.2 mm
046.799		Alignment Pin Ø 2.2 mm
046.800		Depth Gauge Ø 2.8 mm
046.802		Depth Gauge Ø 3.5 mm
046.804		Depth Gauge Ø 4.2 mm

046.460		Implant Adapter for Ratchet, extra short
046.461		Implant Adapter for Ratchet, short
046.470		Implant Adapter for Handpiece, extra short
046.471		Implant Adapter for Handpiece, short

B Module Tissue Level Specific Tools

Add for Standard Plus and Standard Plus Short Implant line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
044.086		SP Profile Drill, short, Ø 3.3 mm
044.088		SP Profile Drill, short Ø 4.1 mm
044.084		SP Profile Drill, short, Ø 4.8 mm
044.575		S/SP Tap, short, Ø 3.3 mm
044.577		S/SP Tap, short, Ø 4.1 mm
044.579		S/SP Tap, short, Ø 4.8 mm

B Module Bone Level Specific Tools


Add for Bone Level Implant Line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
026.2303		BL/TE/NNC Profile Drill, short, Ø 3.3 mm
026.4303		BL Profile Drill, short, Ø 4.1 mm
026.6303		BL Profile Drill, short, Ø 4.8 mm
026.2310		BL/TE/NNC Tap, short, Ø 3.3 mm
026.4310		BL Tap, short, Ø 4.1 mm
026.6310		BL Tap, short, Ø 4.8 mm



A Module Set-up: Standard Plus, Standard Plus Short, Bone Level

A Module
041.761





B Module Set-up: Standard Plus, Standard Plus Short, Bone Level

Parallel Walled Tray 041.779




Implant Adapter for Ratchet
046.460
046.461



Implant Adapter for Handpiece
046.470
046.471




Drill Extender
040.563



Bone Level Tools



Tissue Level Tools



Round Bur
044.004
044.003



Needle Drill
026.0056



Pilot Drill 1
044.210
044.211



Pilot Drill 2
044.214
044.215



Alignment Pin
046.799



Depth Gauge
046.800



Depth Gauge
046.802



Depth Gauge
046.804



Twist Drill PRO
044.250
044.251



Twist Drill PRO
044.254
044.255

TE Order List: Tapered Effect


A Module – Order list

TE Tapered Effect		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device
046.064		Holding Key
041.767		Released Aid Tray
026.2558		Released Aid N for Loxim
026.4558		Released Aid R/W for Loxim
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS screwdriver for Ratchet, extra-short
046.401		SCS screwdriver for Ratchet, short
046.402		SCS screwdriver for Ratchet, long
046.410		SCS screwdriver for Handpiece, extra short
046.411		SCS screwdriver for Handpiece, short
046.412		SCS screwdriver for Handpiece, long
041.762		Grommet Tray, 6 small
044.304		Bone Profiler RN
044.305		Bone Profiler WN
041.763		Grommet Tray, 6 large
80103		48 hour Explantation Device RN/WN
80087		48 hour Explantation Device RN/WN

For additional tools (planning aid, tweezers, AS Screwdriver), please refer to page 5 *A Module full list*.

TE Order List: Tapered Effect

B Module – Order list


TE Tapered Effect		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.778		B Module, Basic Tray
040.563		Drill extender
044.003		Round bur, Ø 2.3 mm
044.004		Round bur, Ø 3.1 mm
026.0056		Needle drill long, Ø 1.6 mm
044.210		Pilot Drill 1, short, Ø 2.2 mm
044.214		Pilot Drill 2, short, Ø 2.8 mm
044.250		Twist Drill PRO, short, Ø 3.5 mm
044.254		Twist Drill PRO, short, Ø 4.2 mm
044.211		Pilot Drill 1, long, Ø 2.2 mm
044.215		Pilot Drill 2, long, Ø 2.8 mm
044.251		Twist Drill PRO, long, Ø 3.5 mm
044.255		Twist Drill PRO, long, Ø 4.2 mm
044.701		TE Profile Drill, short, Ø 3.3 mm
044.705		TE Profile Drill, short, Ø 4.1 mm
044.703		TE Profile Drill, short, Ø 4.8 mm

026.2310		BL/TE Tap, short, Ø 3.3 mm
026.4310		BL/TE Tap, short, Ø 4.1 mm
026.6310		BL/TE Tap, short, Ø 4.8 mm
046.799		Alignment Pin, Ø 2.2 mm, BLT
046.800		Depth Gauge, Ø 2.8 mm, BLT
046.802		Depth Gauge, Ø 3.5 mm, BLT
046.804		Depth Gauge, Ø 4.2 mm, BLT
046.460		Implant Adapter for Ratchet, extra short
046.461		Implant Adapter for Ratchet, short
046.470		Implant Adapter for, Handpiece, extra short
046.471		Implant Adapter for, Handpiece, short




A Module Set-up: Tapered Effect

A Module
041.761



TE B Module Set-up: Tapered Effect

Basic Tray 041.778





Order List: Straumann® PURE Ceramic Implant System

A Module – Order list

Straumann® PURE Ceramic Implant System		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device
046.064		Holding Key
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS screwdriver for Ratchet, extra-short
046.401		SCS screwdriver for Ratchet, short
046.402		SCS screwdriver for Ratchet, long
046.410		SCS screwdriver for Handpiece, extra short
046.411		SCS screwdriver for Handpiece, short
046.412		SCS screwdriver for Handpiece, long
041.763		Grommet Tray, 6 large
031.080		48 hour Explantation Device, Straumann® PURE Ceramic Implant (Monotype) RD*
031.081		48 hour Explantation Device, Straumann® PURE Ceramic Implant (Monotype) ND*
032.053		48h Explantation Device, Straumann® PURE Ceramic Implant RD*

* Note: PURE 48h explantation device are for single use only.

For additional tools (planning aid, tweezers, AS Screwdriver), please refer to page 5 *A Module full list*.



Order List: Straumann® PURE Ceramic Implant System

B Module – Order list

Straumann® PURE Ceramic Implant System		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.778		B Module, Basic Tray
040.563		Drill extender
044.003		Round bur, Ø 2.3 mm
044.004		Round bur, Ø 3.1 mm
026.0056		Needle drill long, Ø 1.6 mm
044.210		Pilot Drill 1, short, Ø 2.2 mm
044.214		Pilot Drill 2, short, Ø 2.8 mm
044.250		Twist Drill PRO, short, Ø 3.5 mm
044.254		Twist Drill PRO, short, Ø 4.2 mm
044.211		Pilot Drill 1, long, Ø 2.2 mm
044.215		Pilot Drill 2, long, Ø 2.8 mm
044.251		Twist Drill PRO, long, Ø 3.5 mm
044.255		Twist Drill PRO, long, Ø 4.2 mm
026.2303		BL/NNC Profile Drill, short, Ø 3.3 mm
026.4303		BL Profile Drill, short, Ø 4.1 mm
026.6303		BL Profile Drill, short, Ø 4.8 mm
026.2310		BL/TE/NNC Tap, short, Ø 3.3 mm
026.4310		BL/TE Tap, short, Ø 4.1 mm
026.6310		BL/TE Tap, short, Ø 4.8 mm

046.799		Alignment Pin, Ø 2.2 mm, BLT
046.800		Depth Gauge Ø 2.8 mm, BLT
046.802		Depth Gauge, Ø 3.5 mm, BLT
046.804		Depth Gauge, Ø 4.2 mm, BLT
046.460		Implant Adapter for Ratchet, extra short
046.461		Implant Adapter for Ratchet, short
046.470		Implant Adapter for Handpiece, extra short
046.471		Implant Adapter for Handpiece, short


C Module – Order list

Straumann® PURE Ceramic Implant System		
Art. No.	Image	Product
041.773		Straumann® Modular, Cassette, C Module, Spare Tools
031.102		Position Indicator, PURE ND, a.h. 4 mm, Ø 2.2 mm
031.103		Position Indicator, PURE ND, a.h. 4 mm, Ø 2.8 mm
031.112		Position Indicator, PURE ND, a.h. 5.5 mm, Ø 2.2 mm
031.113		Position Indicator, PURE ND, a.h. 5.5 mm, Ø 2.8 mm
031.123		Position Indicator, PURE RD, a.h. 4 mm, Ø 2.2 mm
031.125		Position Indicator, PURE RD, a.h. 4 mm, Ø 3.5 mm
031.143		Position Indicator, PURE RD, a.h. 5.5 mm, Ø 2.2 mm
031.145		Position Indicator, PURE RD, a.h. 5.5 mm, Ø 3.5 mm



A Module Set-up: Straumann® PURE Ceramic Implant System


A Module
041.761





B Module Set-up: Straumann® PURE Ceramic Implant System

Basic Tray 041.778






C Module Set-up:

Straumann® PURE Ceramic Implant System

C Module Spare Tools
041.773



Position Indicator	031.103	031.125
031.102	031.113	031.145
031.112		
031.123		
031.143		



BLT Order List: Bone Level Tapered


A Module – Order list

BLT Bone Level Tapered		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device
046.064		Holding Key
041.767		Released Aid Tray
026.2558		Released Aid N for Loxim
026.4558		Released Aid R/W for Loxim
041.767		Release Aid Tray*
026.0073		Release Aid S for Loxim*
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS screwdriver for Ratchet, extra-short
046.401		SCS screwdriver for Ratchet, short
046.402		SCS screwdriver for Ratchet, long
046.410		SCS screwdriver for Handpiece, extra short
046.411		SCS screwdriver for Handpiece, short
046.412		SCS screwdriver for Handpiece, long

041.762		Grommet Tray, 6 small
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3

A Module BLT Ø 2.9 Specific Tools

Add for BLT Ø 2.9 implant line




BLT 2.9 Bone Level Tapered Ø 2,9		
Art. No.	Image	Product
041.763		Grommet Tray, 6 large
026.2048		48 hour tool NC, BLT
026.4048		48 hour tool RC, BLT
041.764		Grommet Tray, 3 small + 3 large*
026.0070		Bone Profiler, oval Platform, BLT SC*
026.0072		48 hour tool SC, BLT SC*
025.0042		SC Implant Adapter for Handpiece, long, BLT SC*
025.0043		SC Implant Adapter for Ratchet, long, BLT SC*


*For Bone Level Tapered Small Crossfit implant line only

For additional tools (planning aid, tweezers, AS Screwdriver), please refer to page 5 *A Module full list*.

BLT Order List: Bone Level Tapered

B Module – Order list

BLT Bone Level Tapered		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.778		B Module, Basic Tray
040.563		Drill extender
044.003		Round bur, Ø 2.3 mm
044.004		Round bur, Ø 3.1 mm
026.0056		Needle drill, long, Ø 1.6 mm
026.0001		BLT Drill short, Ø 2.2 mm
026.2200		BLT Drill, short, Ø 2.8 mm
026.4200		BLT Drill, short, Ø 3.5 mm
026.6200		BLT Drill, short, Ø 4.2 mm
026.0002		BLT Drill, long, Ø 2.2 mm
026.2201		BLT Drill, long, Ø 2.8 mm
026.4201		BLT Drill, long, Ø 3.5 mm
026.6201		BLT Drill, long, Ø 4.2 mm
026.0061		BLT Profile Drill, short, Ø 2.9 mm
026.0003		BLT Profile Drill, short, Ø 3.3 mm
026.0005		BLT Profile Drill, short, Ø 4.1 mm
026.0007		BLT Profile Drill, short, Ø 4.8 mm

026.0063		BLT Tap, long, Ø 2.9 mm, BLT
026.0009		BLT Tap, short, Ø 3.3 mm, BLT
026.0010		BLT Tap, short, Ø 4.1 mm, BLT
026.0011		BLT Tap, short, Ø 4.8 mm, BLT
046.799		Alignment Pin, Ø 2.2 mm, BLT
046.800		Depth Gauge, Ø 2.8 mm, BLT
046.802		Depth Gauge, Ø 3.5 mm, BLT
046.804		Depth Gauge, Ø 4.2 mm, BLT
046.460		Implant Adapter for Ratchet, extra short
046.461		Implant Adapter for Ratchet, short
046.470		Implant Adapter for Handpiece, extra short
046.471		Implant Adapter for Handpiece, short

BLT Order List: Bone Level Tapered

A Module 041.761

Ratchet Tray 041.766

Ratchet
046.119

Torque Control
046.049


Holding Key
046.064

Release Aid Tray 041.767

Release Aid N
026.2558

Release Aid R/W
026.4558

Release Aid S*
026.0073



SC Implant Adapter*
for Ratchet
025.0043

48 hour Tool SC*
026.0072

SCS Screwdriver
for ratchet
046.400
046.401
046.402

SCS Screwdriver
for Handpiece
046.410
046.411
046.412

SCS Screwdriver Set

Grommet Tray 3 small + 3 large
041.764



BL Bone Profiler
026.0022
026.0023
026.0024

Bone Profiler Set

Grommet Tray 6 small
041.762



48 hour Tool
NC/RC
026.2048
026.4048

48h tool Set

Grommet Tray 6 large
041.763



Bone Profiler SC*
026.0070




SC Implant Adapter*
for Handpiece
025.0042

SC Auxiliary Set*

Grommet Tray 3 small + 3 large
041.764

BLT Order List: Bone Level Tapered

Basic Tray 041.778



BLX Order List: Bone Level X


A Module – Order list

BLX		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
066.1100		BLX Torque Control Device
046.064		Holding Key
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS screwdriver for Ratchet, extra-short
046.401		SCS screwdriver for Ratchet, short
046.402		SCS screwdriver for Ratchet, long
046.410		SCS screwdriver for Handpiece, extra short
046.411		SCS screwdriver for Handpiece, short
046.412		SCS screwdriver for Handpiece, long
041.771		Implant Depth Gauge Tray
066.2000		Implant Depth Gauge
041.764		Grommet Tray, 3 small + 3 large
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3

For additional tools, tweezers, AS Screwdriver, please refer to page 5 *A Module full list*.

BLX Order List: Bone Level X


B Module – Order list

BLX		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.777		B Module, Fully Tapered Tray
040.563		Drill extender
044.003		Round bur, Ø 2.3 mm
044.004		Round bur, Ø 3.1 mm
026.0056		Needle drill, long, Ø 1.6 mm
066.1301		X VeloDrill™, short, Ø 2.2 mm
066.1302		X VeloDrill™, short, Ø 2.8 mm
066.1303		X VeloDrill™, short, Ø 3.2 mm
066.1304		X VeloDrill™, short, Ø 3.5 mm
066.1305		X VeloDrill™, short, Ø 3.7 mm
066.1306		X VeloDrill™, short, Ø 4.2 mm
066.1307		X VeloDrill™, short, Ø 4.7 mm
066.1308		X VeloDrill™, short, Ø 5.2 mm
066.1309		X VeloDrill™, short, Ø 6.2 mm
066.1701		X VeloDrill™, long, Ø 2.2 mm
066.1702		X VeloDrill™, long, Ø 2.8 mm
066.1703		X VeloDrill™, long, Ø 3.2 mm
066.1704		X VeloDrill™, long, Ø 3.5 mm
066.1705		X VeloDrill™, long, Ø 3.7 mm

066.1706		X VeloDrill™, long, Ø 4.2 mm, BLX
066.1707		X VeloDrill™, long, Ø 4.7 mm, BLX
046.799		Alignment Pin, Ø 2.2 mm
046.800		Depth Gauge, Ø 2.8 mm
046.801		Depth Gauge, Ø 3.2 mm
046.802		Depth Gauge, Ø 3.5 mm
046.803		Depth Gauge, Ø 3.7 mm
046.804		Depth Gauge, Ø 4.2 mm
046.805		Depth Gauge, Ø 4.7 mm
046.806		Depth Gauge, Ø 5.2 mm
046.807		Depth Gauge, Ø 6.5 mm
066.4201		Implant driver for Ratchet, short, BLX
066.4202		Implant driver for Ratchet, long, BLX
066.4101		Implant driver for Handpiece, short, BLX
066.4102		Implant driver for Handpiece, long, BLX


BLX A Module Set-up: Bone Level X

A Module
041.761



BLX B Module Set-up: Bone Level X

Fully Tapered Tray 041.777



BLX Order List: Bone Level X Guided

A Module – Order list

BLX Guided		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
066.1100		BLX Torque Control Device
046.064		Holding Key
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS screwdriver for Ratchet, extra-short
046.401		SCS screwdriver for Ratchet, short
046.402		SCS screwdriver for Ratchet, long
046.410		SCS screwdriver for Handpiece, extra short
046.411		SCS screwdriver for Handpiece, short
046.412		SCS screwdriver for Handpiece, long
041.771		Implant Depth Gauge Tray
066.2000		Implant Depth Gauge
041.764		Grommet Tray, 3 small + 3 large
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3
041.768		Stop Key Tray
034.006		Stop Key, guided

For additional tools, tweezers, AS Screwdriver refer to page 5 *A Module full list*.

B Module – Order list

BLX Guided		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.780		B Module, Fully Tapered Guided Tray
040.563		Drill Extender
034.010		Mucosa Punch, Ø 3.4 mm
034.011		Mucosa Punch, Ø 4.0 mm
034.012		Mucosa Punch, Ø 4.7 mm
034.215		Milling Cutter, Ø 2.8 mm
034.415		Milling Cutter, Ø 3.5 mm
034.615		Milling Cutter, Ø 4.2 mm
066.1301		X VeloDrill™, short, Ø 2.2 mm
066.1302		X VeloDrill™, short, Ø 2.8 mm
066.1303		X VeloDrill™, short, Ø 3.2 mm
066.1304		X VeloDrill™, short, Ø 3.5 mm
066.1305		X VeloDrill™, short, Ø 3.7 mm
066.1306		X VeloDrill™, short, Ø 4.2 mm
066.1307		X VeloDrill™, short, Ø 4.7 mm
066.1308		X VeloDrill™, short, Ø 5.2 mm
066.1309		X VeloDrill™, short, Ø 6.2 mm
066.1501		X VeloDrill™, medium, Ø 2.2 mm
066.1502		X VeloDrill™, medium, Ø 2.8 mm
066.1503		X VeloDrill™, medium, Ø 3.2 mm

BLX Order List: Bone Level X Guided


066.1504		X VeloDrill™, medium, Ø 3.5 mm
066.1505		X VeloDrill™, medium, Ø 3.7 mm
066.1506		X VeloDrill™, medium, Ø 4.2 mm
066.1701		X VeloDrill™, long, Ø 2.2 mm
066.1702		X VeloDrill™, long, Ø 2.8 mm
066.1703		X VeloDrill™, long, Ø 3.2 mm
066.1704		X VeloDrill™, long, Ø 3.5 mm
066.1705		X VeloDrill™, long, Ø 3.7 mm
066.1706		X VeloDrill™, long, Ø 4.2 mm
066.1707		X VeloDrill™, long, Ø 4.7 mm
046.799		Alignment Pin, Ø 2.2 mm
046.800		Depth Gauge, Ø 2.8 mm
046.801		Depth Gauge, Ø 3.2 mm
046.802		Depth Gauge, Ø 3.5 mm
046.803		Depth Gauge, Ø 3.7 mm
046.804		Depth Gauge, Ø 4.2 mm
046.805		Depth Gauge, Ø 4.7 mm
046.806		Depth Gauge, Ø 5.2 mm
046.807		Depth Gauge, Ø 6.5 mm
066.4401		Guided implant Driver for Ratchet, screw retained
066.4403		Guided implant Driver for Handpiece
066.4404		Guided implant Driver for Ratchet

C Module – Order list

BLX Guided		
Art. No.	Image	Product
041.772		Straumann® Modular Cassette, C Module, Guided surgery
034.291		Drill handle, Ø 2.2 mm, 1 mm/3 mm
034.292		Drill handle, Ø 2.8 mm, 1 mm/3 mm
034.293		Drill handle, Ø 3.2 mm, 1 mm/3 mm
034.294		Drill handle, Ø 3.5 mm, 1 mm/3 mm
034.295		Drill handle, Ø 3.7 mm, 1 mm/3 mm
034.296		Drill handle, Ø 4.2 mm, 1 mm/3 mm
034.284		Drill for Template Fixation Pin, Ø 1.3 mm
034.282		Template Fixation Pin, Ø 1.3 mm
034.298		Template Fixation Pin, Ø 2.8/2.8 mm
034.285		Template Fixation Pin, Ø 5/2.8 mm
034.286		Template Fixation Pin, Ø 5/3.2 mm
034.287		Template Fixation Pin, Ø 5/3.5 mm
034.288		Template Fixation Pin, Ø 5/3.7 mm
034.289		Template Fixation Pin, Ø 5/4.2 mm


BLX A Module Set-up: Bone Level X Guided

A Module 041.761




BLX B Module Set-up: Bone Level X Guided

Fully Tapered Guided Tray 041.780



BLX C Module Set-up: Bone Level X Guided

C Module Guided Surgery 041.772



Notes

Notes

International Headquarters

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

www.straumann.com

© Institut Straumann AG, 2019. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.