





NEODENT KEY RESTORATIVE OPTIONS

One prosthetic connection for all Grand Morse Implants



Neo Manual Screwdriver		GM Cover Screws		GM Healing Abutments							
	104.058 Short		0 mm117.021 2 mm117.022		Profile	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm
	104.060 Medium				Ø 3.3	106.207	106.208	106.209	106.210	106.211	106.212
	104.059 Long				Ø 4.5	106.213	106.214	106.215	106.216	106.217	106.218
					Profile	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm	6.5 mm
					Ø 5.5	106.223	106.224	106.225	106.226	106.227	
				Customizable	Ø 7.0		106.228	106.229	106.230	106.231	106.232

Grand Morse Implant Impression Copings

GM Exact Open Tray	GM Exact Closed Tray
 108.162 Regular 108.163 Long	 108.160 Regular 108.161 Long

Multiple/Single-unit Cement Retained or Single-unit Screw Retained	Single-unit						Multiple-unit	Single-unit						Multiple-unit		Single/Multiple-unit								
	Screw/Cement Retained						Screw Retained		Cement Retained						Removable	Screw Retained								
GM Titanium Block Ø 11.5mm Ø 15.8 mm 135.252 135.253	GM Titanium Base Ø 3.5 Ø 4.5 Ø 5.5 0.8 mm 1.5 mm 2.5 mm 3.5 mm 4.5 mm 5.5 mm 135.202 135.203 135.204 135.205 135.206 135.207 135.214 135.215 135.216 135.217 135.218 135.219 135.235 135.236 135.237 135.238 135.239 135.240 135.208 135.209 135.210 135.211 135.212 135.213 135.220 135.221 135.222 135.223 135.224 135.225 135.241 135.242 135.243 135.244 135.245 135.246						GM Cocr Abutment Ø 3.5/3.75 Ø 4.0/4.3 Ø 5.0/6.0 118.309 118.310 118.311 The set includes one GM CoCr Base, one Titanium Screw and one GM Implant Analog.		GM Exact Abutment 0.8 mm 1.5 mm 115.237 115.238 2.5 mm 3.5 mm 115.239 115.240 4.5 mm 5.5 mm 115.241 115.242		GM Exact Mini Conical Abutment 17°/30 ° 17° 1.5mm 2.5mm 3.5mm 115.249 115.250 115.251 30° 1.5mm 2.5mm 3.5mm 115.252 115.253 115.254		GM Exact Click Universal Abutment Straight 0.8 mm 1.5 mm 2.5 mm 3.5 mm 4.5 mm 5.5 mm Ø 3.3 Ø 4.5 114.566 114.567 114.568 114.569 114.570 114.571 114.578 114.579 114.580 114.581 114.582 114.583 Straight 0.8 mm 1.5 mm 2.5 mm 3.5 mm 4.5 mm 5.5 mm Ø 3.3 Ø 4.5 114.572 114.573 114.574 114.575 114.576 114.577 114.584 114.585 114.586 114.587 114.588 114.589 17° 1.5 mm 2.5 mm 3.5 mm Ø 3.3 Ø 3.5 114.542 114.543 114.544 114.545 114.546 114.547 114.548 114.549 114.550 114.551 114.552 114.553 30° 1.5 mm 2.5 mm 3.5 mm Ø 3.3 Ø 4.5 114.554 114.555 114.556 114.557 114.558 114.559 114.560 114.561 114.562 114.563 114.564 114.565						GM Equator Attachment 1.5mm 2.5mm 102.143 102.144 3.5mm 4.5mm 102.145 102.146 5.5mm 102.147		GM Mini Conical Abutment 0.8 mm 1.5 mm 115.243 115.244 2.5 mm 3.5 mm 115.245 115.246 4.5 mm 5.5 mm 115.247 115.248		GM Micro Abutment 0.8 mm 1.5 mm 115.255 115.256 2.5 mm 3.5 mm 115.257 115.258 4.5 mm 5.5 mm 115.259 115.260	
	20 N.cm														32 N.cm									



Neo Screwdriver

105.132

Note: For all Copings, Healing Abutments and Cover Screws apply 10 N.cm



Hexagonal
Prosthetic Driver

105.137



Manual screwdriver torque
(connects all drivers)

104.050

GM Height
Measurer

128.028

