


T-Fix



Order Number

68-115-004, 1 pack (10 pieces)

can be stored in Bone-Fixation-Set (68-720-002)



 **USTOMED**[®]
INSTRUMENTE

USTOMED INSTRUMENTE
Ulrich Storz GmbH & Co. KG

Bischof-Sproll-Straße 2
78532 Tuttlingen / Germany

Hersteller von Chirurgie- und
Dentalinstrumenten

Phone: +49 (0) 74 61 / 9 65 85 0
Fax: +49 (0) 74 61 / 9 65 85 65

Manufacturer of Surgical and
Dental Instruments

info@ustomed.de
www.ustomed.de

T-Fix

„(in)sane in the membrane“

www.ustomed.de



T-Fix Membrane Screws

designed by Prof. Dr. med. Dr. med. dent. Daniel Rothamel



Competence in Instruments
MADE IN GERMANY



Prof. Dr. med. Dr. med. dent. Daniel Rothamel (MD DMD DDS OMFS MSc)



Education:

1996 – 2001	Dental School, Heinrich-Heine University, Duesseldorf, Germany
2001 – 2007	Medical School, Heinrich-Heine University, Duesseldorf, Germany
2004	Dr. med. dent., magna cum laude
2008	Dr. med., magna cum laude

Professional Experience:

2002 – 2007	Research Associate, Department of Oral Surgery (Prof. Dr. J. Becker) Heinrich-Heine University Düsseldorf, Germany
Feb. – Mar. 2005	Volunteer at the Clinic of Plastic and Reconstructive Surgery in Sankhu, Kathmandu, Nepal
2007	Specialization in Oral Surgery
June – Nov. 2007	Post-doc Student at the Department of Oral Pathology and Molecular Research, Westmead Hospital, University of Sydney, Sydney, Australia (granted by DAAD)
2008 – 2013	Research Associate, Department of Oral and Maxillofacial Plastic Surgery (Prof. Dr. Dr. J. Zöller), University of Cologne, Cologne, Germany
2009	Habilitation, Thesis: "Reconstruction of defects of the alveolar ridge using artificial and autogenous bone blocks and growth factors"
2013 – 2016	Executive senior physician at the Department of Oral and Maxillofacial Plastic Surgery (Prof. Dr. Dr. J. Zöller), University of Cologne, Cologne, Germany
2015 – 2016	Associate Professor for Oral and Maxillofacial Surgery at the Department of Oral and Maxillofacial Plastic Surgery (Prof. Dr. Dr. J. Zöller), University of Cologne, Cologne, Germany
2016 – 2017	Associate Professor and vice head of department at the Department of Oral and Maxillofacial Plastic Surgery (Prof. Dr. Dr. N. Kübler), University of Düsseldorf, Düsseldorf, Germany
since July 2017	Head of the Department of Oral and Maxillofacial Plastic Surgery of Johanniter Hospital Bethesda, Mönchengladbach, Germany

> 100 Publications (71 pubmed-listed, >140 impact points), multiple national and international presentations, various scientific grants and research awards



Ø 1,5 mm, 4 mm length

conical and aggressive thread

specially for rigid or titanium reinforced membranes or meshes

self-cutting and -drilling for (in)sane and easy handling