


STRAUMANN® **VARIOBASE®** XC

Efficiency and flexibility in one.



STRAUMANN® VARIOBASE® XC EFFICIENCY AND FLEXIBILITY IN ONE

Discover the confidence that comes with a system trusted by over 95,000 dental professionals worldwide. The Straumann® Variobase® XC portfolio for the Straumann iEXCEL™ Performance System empowers both dental technicians and restorative dentists to deliver a tailored solution for each unique case with unrivaled flexibility and precision.



RELIABILITY

Deliver fit, structural integrity and retention you can rely on with Straumann® Variobase® XC and its original Straumann implant-abutment connection. The next generation builds on the clinically proven performance of its predecessor, with enhanced abutment-crown bonding strength that reduces the risk of debonding.



EFFICIENCY

Straumann® Variobase® XC helps you achieve uncompromised esthetics with ease from start to finish. Streamline ordering with the Straumann® TorcFit™ connection and customization with the versatile 7 mm high chimney. Skip sandblasting at the end thanks to the retention of the patent-protected laser surface.



FLEXIBILITY

Whether for screw-retained or cement-retained, digital or conventional workflows, the comprehensive and streamlined Straumann® Variobase® XC portfolio enables you to offer customized restorations with ease, even in challenging cases.



STRAUMANN® VARIOBASE® XC

THE NEXT GENERATION OF VARIOBASE®

1 INTEGRATED CUTTING GROOVES

Customize the 7 mm high chimney quickly and precisely.

2 PATENT-PROTECTED LASER-TREATED SURFACE


Enhances cement retention for improved bonding.

3 ANTI-ROTATION NOTCHES

Prevents rotation and increases surface contact for extra stability.

4 RANGE OF GINGIVA HEIGHTS

Multiple options available to meet your patients' needs.





ANGULATION FLEXIBILITY

Supports up to 25° angulation, expanding treatment possibilities.



VERSATILE APPLICATIONS

Suitable for both single- and multi-unit restorations.


SEAMLESS DIGITAL INTEGRATION

Compatible with your preferred digital workflow systems.



CONSISTENT EMERGENCE PROFILES™

Optimize the soft tissue management process with the uniform emergence profile of the Straumann iXCEL™ bone-level abutments.



ENHANCED BONDING FOR A SURE FIT


MACHINED SURFACE

Straumann® Variobase®




LASER-TREATED SURFACE*

Straumann® Variobase® XC



*Patent-protected

ALL THE BENEFITS OF SANDBLASTING. WITHOUT THE RISK.



NOTE: The shortest length (3.5 mm) was used for this laboratory study.

We roughen the abutment surface with a proprietary laser treatment, eliminating the step of sandblasting and helping you to avoid the known risks of the process:

- Altering tolerances or structural components
- Reducing structural strength
- Potential contamination

Enjoy the peace of mind that comes with Straumann® Variobase® XC, knowing Straumann's Swiss engineering and precision manufacturing is completely intact.


THE FOUNDATION OF GENUINE QUALITY

STRAUMANN®


ORIGINAL

The confidence that comes with peace of mind:


- Design harmony between implant and abutment with very high manufacturing precision
- Excellent fatigue test performance
- Long-term clinical evidence
- Guaranteed long-term availability of prosthetic components
- Exceptional service and support including the Straumann Guarantee



COMPATIBLE IS NOT ORIGINAL. WHICH ONE DO YOU TRUST?



Straumann® CARES® Ti Abutment
on Straumann® RC 4.1 Implant



NobelProcera™ CAD/CAM Ti abutment
on Straumann® RC 4.1 Implant

Atlantis™ CAD/CAM Ti abutment on
Straumann® RC 4.1 Implant

LIFETIME GUARANTEE ON STRAUMANN® VARIOBASE® XC ABUTMENTS


- We are so confident in the performance and reliability of Straumann® Variobase® XC that we offer the Straumann Lifetime Guarantee
- For more details, refer to www.straumann.com/guarantee



STRAUMANN® VARIOBASE® XC

DIGITAL WORKFLOW

NOTE: Centralized and conventional workflows are also available.



STRAUMANN® VARIOBASE® XC PORTFOLIO



FOR SINGLE-UNIT RESTORATIONS

		FOR CROWN		FOR CROWN AS	
					
		Gingiva height (mm)	Article number	Gingiva height (mm)	Article number
Straumann iEXCEL™ BL (BLX/BLC)	RB/WB Ø 3.8	1.5	062.5017	Not available	
		2.5	062.5018		
		3.5	062.5019		
	RB/WB Ø 4.5	1.5	062.5020	1.5	062.5009
		2.5	062.5021	2.5	062.5010
		3.5	062.5022	3.5	062.5011
	WB Ø 5.5	0.75	062.5023	0.75	062.5012
		1.5	062.5024	1.5	062.5013
		2.5	062.5059	2.5	062.5060
	WB Ø 6.5	1.5	062.5061	1.5	062.5063
		2.5	062.5062	2.5	062.5064
Straumann iEXCEL™ TL (TLX/TLC)	NT		062.5025		062.5014
	RT		062.5026		062.5015
	WT		062.5027		062.5016

STRAUMANN® VARIOBASE® XC PORTFOLIO


MULTI-UNIT RESTORATIONS




IMPLANT-LEVEL

		FOR BARS/BRIDGES		FOR BARS/BRIDGES AS	
					
		Gingiva height (mm)	Article number	Gingiva height (mm)	Article number
Straumann iEXCEL™ BL (BLX/BLC)	RB/WB Ø 3.8	1.5	062.5065	Not available	
		2.5	062.5066		
		3.5	062.5067		
	RB/WB Ø 4.5	1.5	062.5068	1.5	062.5073
		2.5	062.5069	2.5	062.5074
		3.5	062.5070	3.5	062.5075
Straumann iEXCEL™ TL (TLX/TLC)	WB Ø 5.5	0.75	062.5071	0.75	062.5076
		1.5	062.5072	1.5	062.5077
	NT		062.5003		062.4992
	RT		062.5004		062.4993
	WT		062.5005		062.4994

STRAUMANN® VARIOBASE® XC PORTFOLIO MULTI-UNIT RESTORATIONS


ABUTMENT-LEVEL



FOR CROWN		FOR BARS/BRIDGES		FOR BARS/BRIDGES AS	
					
SRA Diameter (mm)	Article number	SRA Diameter (mm)	Article number	SRA Diameter (mm)	Article number
4.6	023.0038	4.6	023.0040	4.6	023.0041

*Straumann iEXCEL™ BL (BLX/BLC)

STRAUMANN® VARIOBASE® XC – EFFICIENCY AND FLEXIBILITY IN ONE.



Discover
Straumann® Variobase® XC

- Designed for freedom, engineered for confidence
- Trusted by experienced dental professionals,
perfected for today's challenges



491.048/en/A/00 07/25

International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com

© Institut Straumann AG, 2025. All rights reserved.
Straumann® and/or other trademarks and logos from Straumann® mentioned herein are
the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

